

AQUATEC®

Toilet seat raisers

a range of adaptable toileting solutions



Providing dependable toileting support

Adjustable, reliable, comfortable



Yes, you can.®



Introducing the Aquatec range of raised toilet seats

An adaptable range offering supportive toileting solutions.

Providing a reliable and supportive toileting solution for individuals is very important to aid independence.

The dependable Aquatec range offers a secure and stable environment with fixations that support the placement of the raised toilet seats in situ.

The range offers a variety of possible seat heights to accommodate individual sizes and to ensure optimal comfort.

Featuring models that offer supportive armrests for stability and standing transfers, including the swing-away options for ease of functionality.



Aquatec 90

- > Raises the toilet seat by 100 mm.
- > Attached safely to the toilet with two side clamps.
- > Large hand-screws for quick attachment and detachment.
- > Particularly high load capacity: 225 kg.



Aquatec 900

- > Individual adjustment to the user at any time, no tools required:
 - three different seat heights (6 cm, 10 cm, 15 cm),
 - three different angles, lower at the front.
- > Swing-away, ergonomic armrests for additional support when sitting down and standing up.
- > Available with or without armrests.
- > Seat can be folded up in all setting positions – ease of cleaning.



Aquatec 90000

- > Available in three different seat heights: 2 cm, 6 cm, 10 cm.
- > Can be ordered with or without armrests.
- > Swing-away, ergonomic armrests for additional support when sitting down and standing up.
- > Innovative clip system: The lid and toilet seat can be removed from the base; no tools are required.



Technical data

For more comprehensive pre-sales information about this product range, including the product's user manuals, please see your local Invacare website.



Seat width ▾



Seat depth ▾



Seat height ▾



Total width ▾

	Seat width ▾	Seat depth ▾	Seat height ▾	Total width ▾
Aquatec 90	365 mm	395 mm	105 mm (front) 115 mm (rear)	380 mm (with lid) 375 mm (without lid)
Aquatec 900	360 mm	425 mm	60 / 100 / 150 mm	590 mm (with armrests) 430 mm (without armrests)
Aquatec 90000	365 mm	390 mm	20 / 60 / 100 mm	540 mm (with armrests) 380 mm (without armrests)



Total depth ▾



Total height ▾



Width between armrests ▾



Width of opening ▾

	Total depth ▾	Total height ▾	Width between armrests ▾	Width of opening ▾
Aquatec 90	410 mm	135 mm (with lid) 125 mm (without lid)	-	210 mm
Aquatec 900	490 mm	380 mm (with armrests) 180 mm (without armrests) 200 mm (seat to armrests)	510 mm	225 mm
Aquatec 90000	475 mm	230 / 270 / 310 mm (with armrests) 30 / 70 / 110 mm (without armrests) 200 mm (seat to armrests)	450 mm	205 mm



Depth of opening ▾



Weight ▾



Load capacity ▾

	Depth of opening ▾	Weight ▾	Load capacity ▾
Aquatec 90	270 mm	1.6 kg (with lid) 1.1 kg (without lid)	225 kg
Aquatec 900	285 mm	5.0 kg (with armrests) 2.7 kg (without armrests)	120 kg
Aquatec 90000	275 mm	4.4 kg (SH:100 mm) 4 kg (SH:60 mm) 3.6 kg (SH:20 mm) (with armrests)	150 kg



Frame colour ▾

Tel: +41 61 487 70 70
hqeurope@invacare.com
www.invacare.eu.com

Invacare International GmbH
Benkenstrasse 260
4108 Witterswil
Switzerland

