



**CHASSWHEEL**  
Possible.

**NO SUCH THING  
AS UNBEATABLE  
CHALLENGES**



Trusted Solutions, Passionate People



Yes, you can.®

# FOUR X

10 year  
Frame warranty

## STRONG

The sturdy, yet flexible chassis, efficient shock absorption and advanced four-wheel drive form a solid base for smooth movements even on challenging surfaces and areas.

## FLEXIBLE

The Four X is comprehensively equipped with additional features such as seat tilt, balance adjustment, legrest and armrest adjustments to enable an individual to lead as active a lifestyle as possible. Upgrade to lift or stand support to get an even more tailored experience.

## STURDY

Every Chasswheel is manufactured using quality parts and assembled by highly trained experts. Four X comes with the industry's most comprehensive 10 year warranty.

## GOING PLACES

Chasswheel produce top of the range all terrain power chairs for individuals who choose an active lifestyle. Its innovative design features 4-wheel suspension and flexible chassis that makes it ideal for challenging terrain. With Chasswheel, going off the beaten track provides an exciting, all round driving experience.



# FOUR X URBAN



## COMPACT

Urban is compact in size and easily manoeuvres the city streets, shops or office environment.

## VERSATILE

The four wheel drive packs plenty of power and performance, easily tackling curbs and uneven terrain smoothly and comfortably.

## SMOOTH

Its advanced steering system provides a great deal of agility giving a versatile, smooth drive when manoeuvring in confined spaces.

## BUSINESS OR PLEASURE

Chasswheel is the ideal choice for smooth manoeuvring in any given environment, whether travelling around town or venturing further afield. Its sturdy design will tackle most terrains unlike conventional wheelchairs.



## STANDARD FEATURES

- Patented flexible chassis
- Tilt
- Gravity point adjustment
- LED indicators
- LED driving lights
- LED rear lights
- Headrest
- Seatbelt
- Adjustable armrest
- Adjustable legrest
- Recaro Specialist seats



## TECHNICAL SPECIFICATIONS

	FOUR X URBAN	FOUR X
Width:	620 mm	695 mm
Height (normal seating position w/out headrest):	1000 mm	700-850 mm
Length:	920 mm (1170mm footrest down)	1250 mm (footrest down)
Weight with battery:	135 Kg With lift 150 Kg	140 Kg Four X DL Lift 170 Kg Four X DL SSS 170 Kg
Maximum speed:	10km/h	10 Km/h
Drive distance:	~40km	~40km
Batteries:	2 x 73Ah, gel	2 x 73Ah, gel
Motors:	4 x 250W	4 x 250W
Max user weight:	125 kg	125 Kg (160Kg with HD kit)
Seat width:	400-480 mm	400-480 mm
Seat depth:	390-510 mm	390-510 mm
Backrest height:	500-600 mm	500-600 mm
Ground clearance:	6" tyres 70mm	6" tyres 100mm 8" tyres 120 mm
Slope climbing capacity:	20°	20°
Maximum obstacle climbing ability:	150 mm	150 mm
Seating height without cushions;	6" tyres - 500mm With lift - 590 mm	6" tyres - Four X DL 440 mm Four X Lift - 530 mm Four X SSS - 530 mm 8" tyres - Four X DL 465 mm Four X Lift - 555 mm Four X SSS - 555 mm
Gravity point adjustment:	200 mm	250 mm
Tilt angle:	30°	30°
Footrest angle:	80°	80°
Backrest angle:	84°-125°	84°-125°
Turning radius:	1.15m	1.15m
Electronics:	Dynamic DX2 REM420/550/	Four X, Four X DL, Four X Lift Dynamic DX2 Four X SSS - R-Net



Think\_Mobility



Invacare Ltd - Think Mobility

**Invacare Ltd**  
Unit 4 Pencoed Technology Park  
Bridgend CF35 5AQ

T: 01656 776200 F: 01656 776201

e: uk@invacare.com www.invacare.co.uk



**Yes, you can.®**

Chasswheel brochure - 2015 - Pt no. 1589554

© 2015 Invacare International Sàrl. All rights reserved. All information quoted is believed to correct at time of print. Invacare Ltd distribute this product under agreement from Chasswheel, Finland.