

# Softform<sup>®</sup> Trolley Mattress Range



Pressure Redistributing A&E  
Trolley Mattresses



Yes, you can.

The *Invacare **Softform Trolley mattress range*** has been designed to specifically meet the demands of the modern A&E department. Suitable for those 'At Risk' of developing a pressure ulcer, the high specification pressure redistributing foam mattress offers a comfortable support surface and interim pressure area care for patients waiting for treatment in A&E. The mattresses come in a variation of sizes and options for use on a selection of trolley frames and are suitable for X-ray and theatres.

### Effective Pressure Ulcer Prevention

The foam core has independent geometrically-cut surface cells which reduce shear and friction whilst increasing the overall support surface area for the patient, thereby reducing peak pressures.

### Multi stretch PU cover

The cover fabric provides excellent recovery and durability and helps to reduce shear and friction forces. The base of the mattress is a tough coated polyurethane (PU), which protects the mattress against potential bed frame damage.

## Features and Options



### Available in a range of heights

TM20: 50mm  
TM40: 100mm  
TM60: 120mm

### Buckle and Straps

This free of charge optional extra prevents the mattress from slipping whilst on the trolley.

### Inspection zip

All mattresses have a zip to allow for easy inspection. A high frequency welded flap, concealing the zip, is available on the TM60 only.

### TM60 additional options

The 120mm height, TM60, trolley mattress is available with visco elastic (memory) foam, octagonal corners and/or Hook&Loop fastenings, making it suitable for a range of platforms and requirements.

## Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

	Width	Length	Height	Total product weight	Max. user weight	Max. cleaning temperature	Warranty
<b>TM20</b>	560 - 700 mm	1830 - 1930 mm	50 mm	5.5 kg	17 stone/108 kg	80°	1 year
<b>TM40</b>	560 - 700 mm	1830 - 1930 mm	100 mm	7 kg	17 stone/108 kg	80°	
<b>TM60</b>	560 - 700 mm	1830 - 1930 mm	120 mm	7 kg	22 stone/140 kg	80°	

Manufactured to European Standards: BS 5852: Crib 5 and BS 7175: Crib 5.

## Colour

Blue



Tel: +44 61 487 70 70  
hqeurope@invacare.com  
www.invacare.eu.com

Invacare International GmbH  
Neuhofweg 51  
4147 Aesch BL  
Switzerland



Yes, you can.®