

Pronto™ M41



Versatile & compact
mobility for everyday use

LINX®



The *Invacare Pronto M41* is a sleek and highly versatile indoor power wheelchair. It offers an excellent range of seating options to address complex postural needs providing support and comfort. A compact mid wheel drive base makes it easy to manoeuvre in most environments and simple to transport. It is designed to perfectly suit an independent lifestyle – wherever you choose to be.

On the spot

The M41 features a mid-wheel-drive base with very small dimensions. Its compact size ensures a tight turning radius and versatile manoeuvrability in constricted spaces.



Everyday obstacles

Invacare's simplified *SureStep*® technology allows individuals to confidently and comfortably drive over everyday thresholds, keeping six wheels on the ground for maximum stability.

The personal fit

The *Pronto M41* offers many seating options with Invacare's highly adjustable modular concept. *Modulite*™ seating includes tension adjustable upholstery styles as well as Flex3 back options featuring a selection of accessories. The seating system interfaces with *Matrix*® backrest options for optimised support and positioning where needed.



LiNX Smart Technology: Redefining Mobility ▼

Invacare LiNX is our **insight inspired** control system with **advanced technology** that provides a **superb driving experience** for users and allows professionals to **configure and tailor** powerchairs quickly and more intuitively. The smart control system enables programming and maintenance updates to be made **wirelessly** and **in real time**, plus it has **straightforward access** to statistical and diagnostic information, giving complete visibility to professionals.



Touchscreen technology

The REM400 remote has an impressive **3.5" colour touchscreen** that requires **minimal force** to operate. Bringing **modern technology** to the powerchair market, this innovative **interactive remote** can be configured to a users' physical and cognitive ability.



Display only version

The REM500 is a **display only** version of the REM400, retaining the 3.5" interactive touch screen, and allows specialist controls to be used to operate the powerchair.



Built in Bluetooth

The REM400 remote allows users to be connected to multiple devices simultaneously and has built in mouse mover and switch control functionality to be able to control PC, Laptop, smart phone etc.



Features and options ▼



▲ Simplified SureStep

Keeps six wheels on the ground for smooth driving and optimum stability.



▲ Tie-downs

4 integrated, crash-tested tie-downs.



▲ Service position

Tool-free access to critical components.



▲ Recline

30° of recline with a high pivot point for a natural movement.



▲ LN_x

Powered centre mount allows for an active seating position.



▲ Flip-up armrest

The **Modulite™** seat can be combined with a multi-adjustable flip-up armrest for personalised comfort and easier transfers.



Invacare **Modulite™**



▲ Legrest options

Choice of manual settable or powered elevating versions.



▲ Matrix backs

Different types of **Matrix®** backs can be combined with **Modulite** to offer optimal postural support and positioning.



▲ Flex3 back

The Flex3 back support adjusts from 38 cm to 53 cm.



▲ Relaxing position

The **Modulite** seating can be adjusted to a full relaxing position with 30° of tilt and 30° of recline.



▲ LiNX G-Trac® technology

G-Trac ensures the directions indicated by the control are followed by the chair.



◀ Lifter

The seat lifter has a range of 30 cm.



◀ Tilt in space

30° of tilt with centre of gravity shift for those requiring enhanced posture support.

Technical data ▼

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

Seat width ▼	Seat depth ▼	Seat height ▼	Backrest height ▼	Armrest height ▼	Length of footrest ▼	Seat angle ▼

Pronto M41	380 - 530 mm ¹	410 - 510 mm ² 510 - 560 mm ²	420 - 720 mm ³	480 - 660 mm ⁴	245 - 310 mm ⁵ 230 - 360 mm ⁶	290 - 460 mm ⁴	0 - 30°
-------------------	---------------------------	--	---------------------------	---------------------------	--	---------------------------	---------

Backrest angle ▼	Total height ▼	Driving unit ▼	Total length inc. Footrests ▼	Total length without footrests ▼	Max. User weight ▼	Battery capacity ▼

Pronto M41	90 - 120°	1030 - 1090 mm ⁴	610 - 690 mm ⁷	1030 mm ⁴	830 mm	136 kg	2 x 12V 32 Ah
-------------------	-----------	-----------------------------	---------------------------	----------------------	--------	--------	---------------

Turning diameter ▼	Max. Climbable obstacle height ▼	Max. Safe slope ▼	Motor capacity ▼	Driving range ▼	Speed ▼	Total weight ▼

Pronto M41	1075 mm	50 mm	7° / 12.3%	2 x 250 W	up to 21 km ⁸	6 km/h	112 to 128 kg ⁴
-------------------	---------	-------	------------	-----------	--------------------------	--------	----------------------------




Invacare® **Modulite™**



1. Settable in 4 steps with additional 50 mm by T-arm adjustment or -50 mm with in-rigged armrests
2. Infinite adjustment
3. With lifter, +30 mm for sling seat option
4. Depending on configuration
5. T-armrest

6. Flip-up or following armrest
7. Wheelbase width 610 mm
8. According to ISO 7176-4:2008
9. Invacare Modulite™ seat components are cross compatible across the Invacare power wheelchair range and partly with the Invacare Rea® family.

Frame colours ▼



Electric Red



Scan here to access the website.

Tel: +41 61 487 70 70
hqeurope@invacare.com
www.invacare.eu.com

Invacare International GmbH
Benkenstrasse 260
4108 Witterswil
Switzerland



Yes, you can.®