



AQUATEC®

Transfer aids





Aquatec *Trans*

- Extra large seat surface
- Safe: turntable is only unlocked by the weight of the user
- Seat surface turns and slides easily
- Soft, skin-friendly cover: easy-care, wipes clean



Aquatec *Disk*

- Extra large seat surface
- Seat surface turns very easily
- Soft, skin-friendly cover: easy-care, wipes clean
- Within seconds, the **Disk**, the shower stool **Dot** and an adapter are combined to a shower stool with a swivel seat



Aquatec *Step*

- Stable frame
- Feet with non-slip rubber caps
- Non-slip, removable cover, machine-washable

Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

| Order numbers | Aquatec <i>Trans</i> | Aquatec <i>Disk</i> | Aquatec <i>Step</i> |
|--------------------------------|----------------------|---------------------|---------------------|
| Turning and transfer aid, blue | 4.03.002 | | |
| Turning and transfer aid, grey | 4.03.008 | | |
| Rotating disk, blue | | 4.02.002 | |
| Rotating disk, grey | | 4.02.008 | |
| Step, blue | | | 3.00.002 |
| Step, without cover | | | 1584940 |
| Dimensions | | | |
| Diameter | 365 mm | 365 mm | - |
| Length | 516 mm | - | 410 mm |
| Width | - | - | 310 mm |
| Height | 25 mm | 25 mm | 165 mm |
| More info | | | |
| Weight | 2.2 kg | 1.6 kg | 1.5 kg |
| Load capacity | 140 kg | 140 kg | 150 kg |
| Colours | Blue, grey | Blue, grey | Blue, white |

QR
Code

Tel: +41 61 487 70 70
hqeuropa@invacare.com
www.invacare.eu.com

Invacare International GmbH
 Benkenstrasse 260
 4108 Witterswil
 Switzerland

