



INVACARE® TOP END®  
**FORCE™ G**  
HANDCYCLE



The **Invacare® Top End® Force™ G Handcycle** design offers an adjustable, trunk-powered, erect position for low paras, amputees and competitors who can “gut” it out using abdominal, arm, shoulder and chest muscles. The Top End Force G Handcycle features an internal rib center frame design to make it super-stiff. The adjustable back angle (55°-90°), back height, leg rests, crank height and fore/aft seat position allows you to dial in your ultimate power position. Are you ready to force the pace?

## Features

Custom skeletal reinforced aluminium frame  
Adjustable seat length  
Adjustable back angle  
Aluminium oversized rake fork  
Adjustable crank housing  
Multi-position adjustable footrests  
Safety straps  
Travel ready front fork  
Top End V or S crankset

Oval aluminium ergonomic handpedals  
Shimano® 27-speed drive train  
Shimano RAPIDFIRE hands-on shifters/brake  
11-34 SRAM® cassette  
Cable driven upper derailleur shifter  
Upper chainrings: 48-39-24  
Steering dampener  
26" high performance wheels  
High pressure clincher tires

Threaded axles  
Drafting bumper  
Parking brake  
Seat and back cushions  
Chain guard  
Safety flag  
250 lb. weight capacity

## Model #: FRCG



1. Adjustable back angle 55°-90°
2. Drafting bumper
3. Cable operated crank mounted shifter for upper chainrings
4. Shimano® RAPIDFIRE shifter/brake mounted on a handpedal with optional power plate
5. Fork with diagonal and adjustable height and chain guard
6. Travel ready front fork



**Invacare Corporation**  
[www.invacare.com](http://www.invacare.com)

### USA

One Invacare Way  
Elyria, Ohio  
44035-4190  
(800) 333-6900

### Canada

570 Matheson Blvd. E.,  
Unit 8  
Mississauga, Ontario  
L4Z 4G4 Canada  
(800) 668-5324

### Top End

[www.topendwheelchair.com](http://www.topendwheelchair.com)  
Sports and Recreation  
Products  
4501 63rd Circle North  
Pinellas Park, Florida  
33781-5914  
(727) 522-8677  
(800) 532-8677  
Fax: (727) 522-1007

©2014

Invacare Corporation.  
All rights reserved.  
Trademarks are identified by the symbols™ and®.  
All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Shimano and Rapid Fire are registered trademarks of Shimano Inc. SRM is a registered trademark of SRM Service Center Inc. USA. SRAM is a registered trademark of SRAM Corp. Continental is a registered trademark of Continental Tire Corp. Specifications are subject to change without notification.  
Form No. 12-450  
Rev. 12/15  
141160

## Options

Carbon fiber wheels  
SRM® training system  
SRAM grip or trigger shifters  
SRAM cassette choices  
11-23 or 11-26  
53-40-28 Chainring  
Continental® Grand Prix clincher tires  
Continental Sprinter tubular tires  
Heart rate monitor  
Alignment gauge  
Backpack hydration system

Disk brake  
Large or small hand pedals  
Power plates  
Water bottle and cage  
Cordless computer  
Helmet  
Safety lights, mirror  
Additional safety flag  
Click straps  
Crutch holder  
Bike rack

Indoor training roller  
Extra cable kit  
Leg guard attachment  
Extra crank arms  
Tool and tire repair kit  
Wheel bag  
Top End bike jersey

