

Colibri®



reddot award 2014
winner

The colourful route to
effortless, safe drive



The new *Invacare Colibri* is a stylish, simple and colourful micro lightweight scooter, designed for those who enjoy a totally independent lifestyle.

Thanks to the unique patented *Invacare LiteLock™* system the scooter can easily be taken apart without the need for tools and fits neatly into any car boot. Plus, with its small footprint the *Invacare Colibri* offers greater manoeuvrability in tight, or limited spaces. Additionally, you can choose from a range of trendy colour panels to match your style, swapping them according to your needs.



Features and Options



Effortless dismantling

Don't think - do it! Less than a minute is all that's needed to disassemble the *Colibri* for transportation.



Option of a bigger wheel set and battery pack

Perfect for travelling longer distances, on bumpier terrain, such as cobbled streets. Choose where to charge your battery: directly on the *Colibri* or detached from it.



Increased comfort

The soft, padded seat provides comfort during driving. For increased comfort there is also the optional 'easy fit' seat suspension kit.



Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local *Invacare* website.

	Backrest width	Seat depth	Seat to floor height	Backrest height	Driving unit width	Max. Climbable obstacle height	Max. Safe slope	Motor capacity
<i>Invacare Colibri</i>								
	465 mm	400 mm	370 - 420 mm	385 mm	505 mm	45 mm	6° / 10.5%	200 W
	Max. Driving distance	Max. Driving speed	Total length	Empty weight	Max. User weight	Battery capacity	Turning diameter	
<i>Invacare Colibri</i>								
4 wheels	Up to 16 km	8 km/h	1010 mm	44 kg (12 Ah) 49 kg (18 Ah)	136 kg	2 x 12 Ah 2 x 18 Ah	2200 mm	
3 wheels	Up to 16 km	8 km/h	1010 mm	42 kg (12 Ah) 46 kg (18 Ah)	136 kg	2 x 12 Ah 2 x 18 Ah	1940 mm	

Shroud colours



LiteLock™ system (patent pending)



QR
Code

Tel: +41 61 487 70 70
hqeurope@invacare.com
www.invacare.eu.com

Invacare International GmbH
Benkenstrasse 260
4108 Witterswil
Switzerland



© 2018 Invacare International GmbH. All rights reserved. All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation. COLIBRI - EU - 10/2018