



INVACARE® TOP END®
SCHULTE 7000 SERIES
BASKETBALL WHEELCHAIR



The **Invacare® Top End® Schulte 7000 Series Basketball Wheelchair** is faster, lighter, stronger and better than ever. In 2011 it hit the market with new material, an enhanced headtube design and an axel securing system, reducing overall weight by two pounds and improving performance. Each Schulte 7000 Series Wheelchair is custom built to an individual athlete's specifications and is overseen by Paul Schulte, the design engineer for Top End who is a 3 time Paralympian having been a member of the USA National Men's Wheelchair Basketball Team.

Features

Custom-built fixed frame in 7005 aluminium

Adjustable height fixed swivel anti-tip 58mm fifth wheel

Fixed back with fixed or adjustable height

Adjustable tension seat and back upholstery

Fixed axle position with fixed wheelbase

Fixed camber: 15°, 18° or 20°

High performance Sun Components® wheels:

24", 25", 26" or 700C

Four bed offensive wing

Anodized pushrims

High pressure clincher tires

Precision fork caster system

72mm casters

Adjustable height tubular footrest


Powder coat paint

Rear frame protector

Axle securing system

250 lb. weight capacity

Model #: PS7



1. Double swivel anti-tip/ fifth wheel option
2. Precision fork caster system
3. Optional welded fixed tube sideguard
4. Adjustable center-of-gravity option
5. Adjustable height fixed swivel anti-tip
6. 7005 Aluminium heat treating process



Invacare Corporation

www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

Top End

www.topendwheelchair.com
Sports and Recreation
Products
4501 63rd Circle North
Pinellas Park, Florida
33781-5914
(727) 522-8677
(800) 532-8677
Fax: (727) 522-1007

©2016

Invacare Corporation.

All rights reserved.

Trademarks are identified

by the symbols ™ and ®.

All trademarks are owned

by or licensed to Invacare

Corporation unless

otherwise noted. Spinerig

is a registered trademark

of Spinerig Inc. Sun and

Equalizer are a registered

trademark of Sun

Components. Specifications

are subject to change

without notification.

Form No. 12-480

120902

Options

Adjustable center-of-gravity

Spinerig® SLX or Spox wheels

Sun Components Equalizer® wheels

12 tab pushrims

Tapered seat width

Welded seat rigidizer bars with solid seat

Double swivel anti-tip/fifth wheel

Click straps

Double amputee front-end

Recessed platform footrest

Adjustable angle footrest

Clamp-on fixed aluminium side guards

Fixed side guards with welded tube

Fixed aluminium or carbon fiber side guards

Titanium axles

Titanium pushrims

Front frame protectors

Alignment gauge



Yes, you can.®