

# Spectra® Blitz



  
**JUNIOR**

Go and grow with the  
Invacare Spectra Blitz

  
Yes, you can.®

The *Invacare Spectra Blitz* is a versatile power chair for children. Its compact design provides excellent mobility and driving performance. The seat is highly adaptable and allows precise adjustment for optimal comfort and posture management.

## Features and options



Height adjustable back rest.



Adjustable seat width from 25 to 39 cm.



















Rear wheel suspension for better shock absorption.



## Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

	Seat width	Seat depth	Seat height	Backrest height	Backrest angle	Driving unit width	Total length inc. Footrests	Total height
<i>Spectra Blitz</i>	 250, 300, 350 mm <sup>1</sup>	 250 - 380 mm <sup>2</sup>	 450 mm	 370 - 440 mm	 -5° - 30°	 560 mm	 940 mm	 910 mm

	Turning diameter	Max. Climbable obstacle height	Max. Safe slope	Battery capacity	Motor capacity	Driving range	Speed	Total weight
<i>Spectra Blitz</i>	 1240 mm	 50 mm	 14.1% (8°) <sup>3</sup>	 2 x 40 Ah	 2 x 180 W	 up to 21 km <sup>4</sup>	 6 km/h	 79.5 kg (incl. batteries) <sup>5</sup>

Max. User weight




*Spectra Blitz* 75 kg

1. Additional 40 mm armrest adjustment
2. In 6 steps
3. The maximum safe slope depends on the configuration. For further details please see the user manual.
4. According to ISO 7176-4:2008
5. Depending on configuration



## Frame colour



Tel: +41 61 487 70 70  
[hqeuropa@invacare.com](mailto:hqeuropa@invacare.com)  
[www.invacare.eu.com](http://www.invacare.eu.com)

Invacare International GmbH  
 Benkenstrasse 260  
 4108 Witterswil  
 Switzerland

