

matrix[®]
SEATING SERIES

Flo-tech™ Lite Visco



Low profile visco elastic low memory
foam cushion offering exceptional comfort



Yes, you can.®

The *Invacare **Matrix Flo-tech Lite Visco*** cushion is made of high-density visco elastic foam, allowing it to conform to the user's shape. This affords exceptional levels of comfort and increases pressure redistribution. This low profile cushion format is ideal for both chair and wheelchair use, facilitating side transfers.



Features and Options



Ingress Resistant Cover

Cover is two-way stretch, water-resistant and vapour-permeable. Cover features an anti-slip base with carry handle.



Available in 61 sizes

The *Invacare **Matrix Flo-tech Lite Visco*** is available in a range of 61 sizes, from 250 x 250 to 480 x 530 mm. (See size chart for details).



Low memory foam designed to maximise body contact area

The visco-elastic foam conforms to the user's body, maximising the contact area, delivering exceptional levels of comfort and optimising pressure redistribution.

Three foam densities available

Available in three densities to accommodate a range of user weights – Soft (up to 32 kg), Medium (up to 51 kg) and Firm (up to 100 kg).

Contoured seating surface with pommel and side support

Offers mild lateral support and pelvic stability to enhance the user's posture.

Technical data


For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

Various sizes	Height	Total product weight	Max. user weight ¹	Max. cleaning temperature
A range of 61 sizes	78 mm	1.26 kg	100 kg	80°

Manufactured to European Standards:
Interior Foam: BS EN 1021-1, BS EN 1021-2,
BS 5852: 1982 Crib 5

1. Guide maximum user weight limits are set assuming that the cushion selected is of appropriate width and depth for the intended user/s.

Colours



Tel: +44 61 487 70 70
hqeurope@invacare.com
www.invacare.eu.com

Invacare International GmbH
Neuhofweg 51
4147 Aesch BL
Switzerland



Yes, you can.®