

# Pronto™ Air



A new breed of pure mobility



The *Invacare Pronto Air* Personal Transporter redefines pure power mobility for the consumer. It is a culmination of years of defining Invacare's approach to mobility, with the delivery being innovation and user experience. It's driving as it's meant to be – on the driver's terms.

### Innovation and design

Efficiency of energy and space is at work with Invacare's in-line motor technology emphasizing travel range and allowing for a smaller chassis construction. The Personal Transporter allows the freedom to travel without expensive additional equipment to transport. Most mid-size and larger automobile trunks will easily accommodate the **Pronto Air PT**.

### Features and options



#### Pocket storage bag

The **Pronto Air PT** comes with a bag for personal belongings which secures to the **MyBody** seating shell.



#### Battery pack

The battery pack easy disconnects and reconnects automatically.



#### Transport

Designed with two grip points on the footplate and rear chassis for one or two person lifting.

### Invacare MyBody Seating

Inspiration is everywhere in the innovative design, but it is especially evident in **MyBody** Seating. Invacare designed **MyBody** Seating to promote better posture through contoured, pressure relieving foam and natural points of rest which may result in less fatigue at the end of a fulfilling, accomplished day.

### Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

	Seat width	Seat depth	Seat height	Backrest height	Armrest height	Backrest angle	Driving unit width	Total height
<b>Pronto Air PT</b>	460, 510 <sup>1</sup> mm	460 - 510 mm	535 - 635 <sup>2</sup> mm	535 mm	155 - 230 mm	8° - 32°	600 <sup>3</sup> mm	1070 - 1170 mm
	Total length inc. Footrests	Total weight	Max. User weight	Battery capacity	Turning diameter	Max. Climbable obstacle height	Max. Safe slope	Speed
<b>Pronto Air PT</b>	960 mm	34 <sup>4</sup> / 25 <sup>5</sup> / 13.6 <sup>6</sup> kg	113 kg	2 x 17 Ah	1484 mm	50 mm	3°	6 km/h
	Driving range	<ol style="list-style-type: none"> <li>Additional 76 mm per armrest</li> <li>In 25 mm increments</li> <li>Wheelbase width</li> <li>Wheelbase w/o batteries</li> <li>Seat</li> <li>Battery pack</li> </ol>			Shroud colour			
<b>Pronto Air PT</b>	up to 16 km							

Tel: +41 61 487 70 70  
[hqeuropa@invacare.com](mailto:hqeuropa@invacare.com)  
[www.invacare.eu.com](http://www.invacare.eu.com)

Invacare International GmbH  
 Benkenstrasse 260  
 4108 Witterswil  
 Switzerland

