



INVACARE® TOP END®  
**ELIMINATOR™ OSR**  
RACING WHEELCHAIR



The **Invacare® Top End® Eliminator™ OSR Racing Wheelchair** is designed for speed. A super-stiff, roll-formed horizontal maintube has been combined with a relaxed headtube angle for exceptional stability and aerodynamics, making it the choice of champion racers like Krige Schabort, Kurt Fearnley, Shirley Reilly and Amanda McGrory. So, if you have your sights set on winning, the Eliminator OSR Racing Wheelchair will help you get to the finish line first.

## Features

6061 T6 oval tube custom aluminium frame  
Integral wrap-around fenders  
3 cage designs: Open-V, U or I-cage  
Kneeling, amputee or traditional upholstery  
700C or 26" high performance rear wheels  
Caliper hand brake  
Aluminium side panels  
Straight axle tube with rear wheel alignment  
Aluminium fork  
Precision track compensator system

20" front wheel with tubular tire  
Rubber-coated 3/8" diameter aluminium rims  
Aero bars  
Powder coat paint  
250 lb. weight capacity



1. I-cage
2. Open-V cage shown with cut-out axle option
3. Oversized rake for maximum speed and cornering
4. Precision track compensator system



**Invacare Corporation**  
[www.invacare.com](http://www.invacare.com)

**USA**  
One Invacare Way  
Elyria, Ohio  
44035-4190  
(800) 333-6900

**Canada**  
570 Matheson Blvd. E.,  
Unit 8  
Mississauga, Ontario  
L4Z 4G4 Canada  
(800) 668-5324

**Top End**  
[www.topendwheelchair.com](http://www.topendwheelchair.com)  
Sports and Recreation  
Products  
4501 63rd Circle North  
Pinellas Park, Florida  
33781-5914  
(727) 522-8677  
(800) 532-8677  
Fax: (727) 522-1007

©2014  
Invacare Corporation.  
All rights reserved.  
Trademarks are identified  
by the symbols ™ and  
©. All trademarks are  
owned by or licensed  
to Invacare Corporation  
unless otherwise noted.  
Corima is a registered  
trademark of Corima,  
S.A. Specifications are  
subject to change without  
notification.  
Form No. 01-192  
Rev 12/15  
141155

## Options

Carbon fiber rear wheels  
Corima™ carbon fiber front wheel  
Alignment gauge  
Wheel bag  
Solid seat/kneel plate  
Fixed footrest  
Cushion  
Cordless speedometer  
Helmet

