

Atlas Lite



Simple and effective



The **Atlas Lite** from *Invacare* is a simple wheelchair that's made from strong lightweight steel and capable of supporting loads of up to 120 kg. Added to that, it's easily adjustable to maximise your patient's comfort and easily transportable with quick release rear wheels. What's more, it's available at a competitive price. So if you're looking for a straightforward wheelchair that's practical, economical and convenient, then look no further than the **Atlas Lite**.



Getting around

The overall mobility of a wheelchair is essential, that's why the *Invacare* **Atlas Lite** is easy to maneuver, while being pushed or self-propelled, indoor or outdoor. It even remains mobile in tight spaces or crowded areas because of its small dimensions and compact frame, proving that its size is never an issue.

Easy to transport and store

Being light weight, compact, and with quick release rear wheels, the **Atlas Lite** is very easy to handle, transport and store. Simply take off the rear wheels, pack it up and reassemble it.

Sit in comfort

At Invacare, comfort is imperative in the design of our products. That's why the **Atlas Lite** also comes as a configurable version with a 30° reclining backrest and elevating legrests, so that you can optimise your patient care.

Features and options



Transit version

With 30° reclining backrest and elevating legrests.



Folding backrest

Reduce overall height for storage and transportation.



Hub brakes

Attendant operated, including push to lock brakes.



Armrests

① Sport, ② desk fixed or height adjustable and ③ full length fixed or height adjustable.

Features and options



Fully equipped version

Includes drip stand, anti-tippers and cane holder.



Rear wheels

① 24" spoked wheels (pneumatic or solid tyre), ② 24" plastic wheels (pneumatic or solid tyre) and ③ 12" wheels (pneumatic or solid tyre).



Headrest extension



Dual handrim version

Ideal for hemi users. Can be mounted on right or left side.



Headrest

Fixed on the backbrace and adjustable in height.



Toilet seat



Plastic table tray























Elevating legrest

Adjustable in angle from 0° to 90°.






Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

| | | | | | | |
|--|--|--|--|---|---|---|
| Seat width | Seat depth | Seat angle | Seat height | Backrest angle | Backrest height | Front frame angle |
|  |  |  |  |  |  |  |
| Invacare Atlas Lite 380 / 410 / 430 / 450 / 480 mm | 420 mm | 3° | 485 / 510 mm | 8° / 10° / 0° - 30° | 420 mm | 80° / 0° - 90° |
| Length of footrests | Armrest height | Driving unit width | Total width, folded | Total length inc. Footrests | Total length without footrests | Total height |
|  |  |  |  |  |  |  |
| Invacare Atlas Lite 420 - 500 mm | 220 mm | Seat width + 205 / 235 mm ¹ | 290 / 320 mm ¹ | 1150 mm | 780 mm | 920 - 965 mm |
| Max. Safe slope | Radius | Max. User weight | Total weight | Mass of the heaviest part |  | |
|  |  |  |  |  | | |
| Invacare Atlas Lite 18,5° | 870 mm | 120 kg | 15,5 kg | 9 kg | | |

1. With hub brakes


Frame colours

| | | |
|---|---|---|
| Chrome | Silver Grey (RAL 9006) | Black (RAL 9005) |
|  |  |  |
| Bordeaux (RAL 3005) | Green (RAL 6004) | |
|  |  | |

Please notice colors may vary slightly from those displayed above.

Upholstery

Black Nylon



Scan here to access the website.

Tel: +41 61 487 70 70
hqeurope@invacare.com
www.invacare.eu.com

Invacare International GmbH
Benkenstrasse 260
4108 Witterswil
Switzerland



© 2018 Invacare International GmbH. All rights reserved. All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation. Atlas Lite - EU - 10/2018

