

Ben9_{NG}



Simple, practical and modular



Yes, you can.

The *Invacare Ben 9_{NG}* is the entry level chair in the Invacare manual range. This next generation of the well known *Ben 9 Plus* and *Duo* combines a lightweight construction with a modular concept that will allow you to change the propulsion type within a few minutes. Transit and self-propulsion versions are available both in standard or heavy duty versions for users weighing up to 160 kg.

Modularity

Invacare has designed the *Ben 9_{NG}* to make life easier. As well as all the well known key features that have given the *Ben 9* such a good reputation, the next generation proposes a new rear wheel support that allows you to change the propulsion type of your chair from a self-propulsion to a transit, or vice versa. This increases adaptability to meet your individual needs.



Added value

The *Ben 9_{NG}*'s standard features include a swing away legrest and detachable armrests for easier transfer. The lightweight folding frame and backrest reduce the overall size when folded for transportation or storage. A wide range of options is available, from contoured cushions to extended backrests and headrests for added support and comfort. Simply choose the options that suit your individual needs and tastes.

Features and Options



Modular rear wheel position

Thanks to the removable rear wheel support you can transform your *Ben 9_{NG}* self-propulsion into a transit version and vice versa.

Compactness

The lightweight frame, backrest and legrests fold away easily to make the chair very compact for storage or transportation.

Heavy Duty version

A heavy duty version is available, ideal for users weighing up to 160 kg user, in both self-propulsion and transit versions. Available in seat widths of 510, 560 and 610 mm with seat depths of 430 or 460 mm.



Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

Invacare <i>Ben9_{NG}</i> 380 - 610 mm	405 / 430 / 455 mm	5°	480 mm	8°	430 mm	220 mm	330 - 430 mm
Invacare <i>Ben9_{NG}</i> Seat width + 200 / 130* mm	1040 / 910 ¹ mm	780 / 670 ¹ mm	965 / 720 (folded) mm	11,5°	965 / 720 (folded) mm	127 / 160 ² kg	19,3 / 17 ¹ / 22 ² kg

Frame colour

Black



Upholstery

Black nylon



Please notice colors may vary slightly from those displayed above.



1. Transit version
2. Heavy Duty version

Tel: +41 61 487 70 70
hqeurope@invacare.com
www.invacare.eu.com

Invacare International GmbH
Neuhofweg 51
4147 Aesch BL
Switzerland



© 2023 Invacare International GmbH. All rights reserved. All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation. EU - 07/23

Yes, you can.®