

# Cosy Up



Great comfort, great design!



Yes, you can.

The **Cosy Up** armchair offers maximum comfort thanks to the cushioned one-piece seat, back and leg cushion. Modern and tasteful, the armchair sits nicely in any home or institutional setting.



#### Depth adjustable backrest

Thanks to the two fixing points, the back of the **Cosy Up** is adjustable in depth, making it suitable for a range of people. By having the depth of the chair set to 450 mm, it makes it suitable for shorter individuals, while the 540 mm deep option is ideal for taller users, providing support where its needed. This adjustable features comes as standard on the armchair.



#### Two seat widths

The **Cosy Up** armchair comes as standard with two different armrest linings, giving the chair two seat width options. Both liners fit easily to the arms with zips and strap fasteners, making them easy to remove or replace.

#### Side pouches


The armchair has a side pouch on each side for things like magazines or newspapers.



The thicker lining gives armchair a seat width of 450 mm.



The thinner lining gives the armchair a seat width of 530 mm.



## Features and Options



#### Removable covers

The armrests, backrest and seat unit covers are attached to the armchair with zips and strap fasteners, making them easy to remove. They can also be machine washed at 30° or 40° depending on the fabric.



#### The control always at hand...

The armrest liners also double up as a place to store and easily access a TV remote.

## Options and Configurations




#### Hand control for single motor

Easy to use ergonomic hand control with two illuminated buttons to operate the armchair.



#### Integrated remote with USB port for single motor

Intuitive and discrete remote, placed to the right. To this integrated remote is added a USB port to recharge a mobile phone or a tablet.



#### Hand control for two motor option

Easy to use ergonomic hand control with four buttons and two specially designed synchronisation buttons.



#### Close to the wall

Each **Cosy Up** with 1 motor can be placed close to a wall without preventing the perfect use of the armchair. Thanks to an optimized conception, developed for small spaces, only 9 cm are needed between the wall and the **Cosy Up** to allow a complete recline position.



## Features and Options



NEW

### Fire retardant base

A new fire retardant base is now available. In compliance with the EN 1021-1 and 1021-2 norms, the PVC fire retardant base can be used with the EasyClean fabric to have a totally fire retardant **Cosy Up**.



NEW

### Padded headrest

For a better comfort during all your rest times. Available in accessory: ref. 1650214



NEW

### New Graphique fabric

To have an armchair always more comfortable, the Graphique set of cushion is available on the **Cosy Up** (ref. 1650213). This design fabric is very easy to clean and fits to every environment.



NEW

### EasyClean Yellow



### EasyClean Blue



### EasyClean Taupe

### EasyClean fabric: clean stains only with water !

Anti-bacterial, EasyClean fabrics are treated to prevent the stains from seeping into the fabric. The EasyClean fabric, allows to clean coffee or chocolate stains only with a clean cloth and water to preserve your armchair.

Convenient and modern, this anti-stain fabric prevent mould and is in compliance with the EN 1021-1 and 1021-2 fireproof norms. Washable in the washer until 40° and conform to the dryer, EasyClean cushions are trendy and comfortable. All these qualities make **Cosy Up**, the perfect armchair for a use in collectivity or at home !



## Features and Options

### Cosy Up: mix, pick and change colour following your mood !

With its 6 colours set of cushion and its 4 bases available, the **Cosy Up** is THE configurable armchair. Change the cushion as often as you change your mind and transform your decoration in a few minutes. **Cosy Up**'s configurability allows to the dealer, in a reduce space, to present many colours changing only the set of cushion.

You can order a complete armchair (base + set of cushion) but also create the **Cosy Up** you want choosing a base and one or many separate set of cushion: now, you can change whenever you want!

## Bases



### 1 Motor PVC Taupe base

Ref. 1650169



### 1 Motor PVC Taupe base, conform with fire proof EN 1021-1 and 1021-2 norms

Ref. 1650207



### 1 Motor PVC Taupe base, with integrated remote and USB port

Ref. 1650206



### 2 Motors PVC Taupe base

Ref. 1650205

## Cushions



NEW

### EasyClean Yellow

Conform with fire proof EN 1021-1 and 1021-2 norms  
Ref. 1650210



NEW

### EasyClean Blue

Conform with fire proof EN 1021-1 and 1021-2 norms  
Ref. 1650211



NEW

### EasyClean Taupe

Conform with fire proof EN 1021-1 and 1021-2 norms  
Ref. 1650212



NEW

### Graphique fabric

Ref. 1650213



### Microfiber Taupe

Ref. 1650208



NEW









### PVC/Polyester Taupe










Ref. 1650209



## Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please refer to your local Invacare website.

	Recline position	Rise position	Lying position	Connection between the seat and legrest	Seat depth	Seat to floor height	Seat width	Max. backrest recline
								
<b>1 Motor Close to the wall</b>	◆	◆		◆	450 ou 540 mm	480 mm	450 ou 530 mm	120°
<b>2 Motors</b>	◆	◆	◆	◆	450 ou 540 mm	480 mm	450 ou 530 mm	135°

	Total height	Total width	Total depth (seated position)	Total depth (recline position)	Wall clearance for maximum backrest recline	Armrest to floor height	Max. user weight	Total weight	Foam density
									
<b>1 Motor Close to the wall</b>	1050 mm	770 <sup>1</sup> / 800 mm	910 mm	1520 mm	90 mm	640 <sup>1</sup> / 700 mm	125 kg	57 kg	Seat 28kg/m <sup>3</sup> Armrests 18kg/m <sup>3</sup>
<b>2 Motors</b>	1050 mm	770 <sup>1</sup> / 800 mm	910 mm	1720 mm	460 mm	640 <sup>1</sup> / 700 mm	125 kg	57 kg	Seat 28kg/m <sup>3</sup> Armrests 18kg/m <sup>3</sup>

- Without armrest liners
- Taupe = Brown-grey

## Complete Cosy Up

**1644970** : 1 Motor - PVC/Polyester base, Taupe colour - Microfiber Taupe cushion  
**1650215** : 1 Motor - PVC/Polyester base, Taupe colour - PVC/Polyester Taupe cushion  
**1650227** : 1 Motor - PVC/Polyester base, Taupe colour - EasyClean Taupe antibacterial cushion  
**1650216** : 1 Motor - PVC/Polyester base, Taupe colour - EasyClean Yellow antibacterial cushion  
**1650217** : 1 Motor - PVC/Polyester base, Taupe colour - EasyClean Blue antibacterial cushion  
**1650228** : 1 Motor - PVC/Polyester base, Taupe colour - Graphique fabric cushion

**1650229** : 1 Motor - PVC/Polyester base with integrated remote and USB port, Taupe colour - Microfiber Taupe cushion  
**1650230** : 1 Motor - PVC/Polyester base with integrated remote and USB port, Taupe colour - PVC/Polyester Taupe cushion  
**1650233** : 1 Motor - PVC/Polyester base with integrated remote and USB port, Taupe colour - Coussin Antibactérien EasyClean Taupe  
**1650231** : 1 Motor - PVC/Polyester base with integrated remote and USB port, Taupe colour - EasyClean Yellow antibacterial cushion  
**1650232** : 1 Motor - PVC/Polyester base with integrated remote and USB port, Taupe colour - EasyClean Blue antibacterial cushion  
**1650234** : 1 Motor - PVC/Polyester base with integrated remote and USB port, Taupe colour - Graphique fabric cushion

**1650235** : 1 Motor - 100% Fireproof - PVC/Polyester base conforms to fireproof norms, Taupe colour - EasyClean Taupe antibacterial cushion  
**1650236** : 1 Motor - 100% Fireproof - PVC/Polyester base conforms to fireproof norms, Taupe colour - EasyClean Yellow antibacterial cushion  
**1650237** : 1 Motor - 100% Fireproof - PVC/Polyester base conforms to fireproof norms, Taupe colour - EasyClean Blue antibacterial cushion

**1644971** : 2 Motors - PVC/Polyester base, Taupe colour - Microfiber Taupe cushion  
**1650238** : 2 Motors - PVC/Polyester base, Taupe colour - PVC/Polyester Taupe cushion  
**1650241** : 2 Motors - PVC/Polyester base, Taupe colour - EasyClean Taupe antibacterial cushion  
**1650239** : 2 Motors - PVC/Polyester base, Taupe colour - EasyClean Yellow antibacterial cushion  
**1650240** : 2 Motors - PVC/Polyester base, Taupe colour - EasyClean Blue antibacterial cushion  
**1650242** : 2 Motors - PVC/Polyester base, Taupe colour - Graphique fabric cushion

Tel: +44 61 487 70 70  
 hqeuropa@invacare.com  
 www.invacare.eu.com

Invacare International GmbH  
 Neuhofweg 51  
 4147 Aesch BL  
 Switzerland

