

matrix[®]
SEATING SERIES

Elan Headrests



Adjustable support with comfort



The elegant *Invacare Matrix Elan* takes our seating range to the **next level** with its **compact, lightweight** and **adjustable** design. The headrest offers **superb support** and **comfort** and compliments a full range of paediatric through to adult applications. Key features include:

- **Highly adjustable hardware**
- **Range of pad options**



Highly adjustable hardware

- The innovative hardware on the Elan headrest has **six points of adjustment** to provide infinite positioning opportunities, including a **swing away function**, to meet user needs
- The headrests can be set up **quickly** and **effectively** with a **single tool** while the **multi-disk locking** feature provides superior holding strength and reliability
- The **tool-less adjustment levers** make it easy to change the position of the headrest
- There are two post options available, the **standard** adult (380 mm) or the **mini** (280 mm), which can be **fitted retroactively** to **Matrix** receivers and are compatible with a number of other headrest receivers on the market



Range of pad options



3 Ellipse pad options

Can be formed to fit the shape of the users head for posterior and lateral support. Available in 3 sizes: 150 x 90 mm, 250 x 130 mm or 360 x 130 mm.



2 Occipital pad options

Provides contoured support to the occipital and sub-occipital areas of the head. Available in two sizes: 230 x 180 mm or 310 x 200 mm.



4-point support

Provides 4 points of support to the occipital and temporal areas of the head. Available in size 280 x 260 mm.

Features

2 tailored cover options

The standard pad cover (Startex fabric) comes with a dual zip and is moisture resistant and breathable. There is also the option of an infection control cover (reverse Startex). Both cover options are easy to remove.

Dual density foam

The high density moulded foam provides **immersion, comfort and protection** for the user and prevents bottoming out. The *Invacare **Matrx*** pads are **45% lighter** than most other pads on the market and surround the edges of the backing plates.



Tool-less adjustments

- 1 Adjustment levers allow easy position changes. Fixed screws provided if set position desired.
- 2 Attachment knob has a slim, ergonomic design to allow easy headrest removal or height adjustment.



Interface insert

An interface insert is available to allow connection to the *Invacare **Flex 3*** headrest receivers.



Multi-disk locking feature

The innovative locking feature provides superior holding strength and reliability.



Adjustable aluminium plates

Can be formed in order to meet the needs of the user on all pad options.



Family of seating products

Invacare can now offer a complete rehabilitation solution: cushion, back support, belt/harness and headrest.

Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

Hardware Options

Elan Standard Hardware	380 mm Vertical Post
Elan Mini Hardware	280 mm Vertical Post


Pad Options

150 mm Ellipse Pad	150 x 90 mm
250 mm Ellipse Pad	250 x 130 mm
360 mm Ellipse Pad	360 x 130 mm
230 mm Occipital Pad	230 x 180 mm
310 mm Occipital Pad	310 x 200 mm
4 Point Pad	280 x 260 mm




Cover options

Startex



Infection control



QR Code



Tel: +41 61 487 70 70
hqeuropa@invacare.com
www.invacare.eu.com

Invacare International GmbH
 Benkenstrasse 260
 4108 Witterswil
 Switzerland

