

**AQUATEC®**

# Ocean Vip / Dual-Vip



Well thought-out details for  
special demands



The tilt-in-space shower commode *Aquatec Ocean Vip* provides plenty of sophisticated details, like the effortless height adjustment and the advanced tilting function. Therefore it assists caregivers in treating severely disabled users who may suffer from spina bifida, multiple sclerosis, para-/tetraplegia in a more gentle and efficient manner. The *Aquatec Ocean Dual-Vip* features full tilt-in-space functions in combination with an adjustable backrest. Care is made easy, even for people with the need of a particular backrest angle.



### Tilting Function

---

The *Aquatec Ocean Vip's* and *Ocean Dual-Vip's* tilting function demands minimum power due to its advanced geometrics. The stepless operation is initiated by an ergonomic lever activating two gas springs to swing the chair backwards (up to 35°). Although both stainless steel gas springs are powerful enough to provide assistance for users up to 150 kg they require only little space within the chair and are placed remotely in order to allow the passing over of large institutional toilets.

### Reclining backrest

---

The reclining backrest of the *Ocean Dual-Vip* is adjustable to backrest angles of 10°- 50°. The backrest can either be retrofitted to existing *Ocean Vip* or be bought as a complete package. Both sets include two anti-tippers for increased safety even in the most extreme positions.

### Adjustment

---

One of the innovative details of the *Ocean Vip* and *Ocean Dual-Vip* is the unique height adjustment. It makes transfer of patients easy as the chair has always the right individually adjusted seating height. The adjustment is a matter of seconds.

### Compatibility

---

For severely disabled persons with reduced lateral stability, standard seating does often not provide enough support. The *Ocean Vip* is compatible with the most common aftermarket accessories and can therefore be equipped with laterals, armsupports or special headrest pads.

## Features and options



### Headrest included

The comfortable headrest can be adjusted steplessly in a wide horizontal and vertical direction.



### Directional lock

One of the four high-quality castors is equipped with a directional lock to help pushing the chair along straight hospital-floors.



### Height adjustable armrests

The armrests can be attached to the backrest in two different heights.



### Adjustable backrest cover

The machine washable backrest cover can be adjusted in a wide range to achieve additional lateral support even without accessories.



### Soft seat with key cut

An even better seating comfort for sensitive patients can be achieved with the soft seat.



### Armrest fixation

The automatically engaging armrest fixation keeps the armrests locked in up and down position.



### Ergonomic lever

The ergonomic lever which operates the tilting function shows advanced design to ease the daily use.



### Reclining backrest

Two stainless steel gas-springs make the step less adjustment of the reclining backrest easy.



### Anti-Tippers

Anti-Tippers are included in **Dual-Vip** sets but can also be bought separately for additional safety.

## Technical data


For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

	Seat height	Seat width	Total width	Total depth	Total height (w/o adjustable headrests)	Tilt angle backrest	Tilt angle seat
<b>Ocean Vip</b>	510 - 610 mm	500 mm	595 mm	1060 mm	1085 - 1185 mm	-	0° - 35°
<b>Ocean Dual-Vip</b>	510 - 610 mm	500 mm	625 mm	1060 mm	1085 - 1185 mm	10° - 50°	0° - 35°

	High load capacity	Product weight
<b>Ocean Vip</b>	150 kg	21.5 kg
<b>Ocean Dual-Vip</b>	150 kg	28.3 kg

## Colour

White



QR Code

Tel: +41 61 487 70 70  
[hqeuropa@invacare.com](mailto:hqeuropa@invacare.com)  
[www.invacare.eu.com](http://www.invacare.eu.com)

Invacare International GmbH  
 Benkenstrasse 260  
 4108 Witterswil  
 Switzerland

