





Adjustments and accessories

Our Ocean Ergo range offers an advanced solution in shower chairs. Each chair is designed to offer high levels of adjustments and configurability as standard. For further customisation, our latest range of accessories provides solutions for extra comfort, support and safety.



Footplates

- Angle adjustable, flip-up and textured footplates with heel loops
- No tools required

Footrests

- Height adjustable, swing away and detachable footrests
- No tools required

Armrests

- Ergonomic, swing away armrests can be adjusted to
- different heights

Seat height

- Seat to ground height can be adjusted between 440-600 mm* or 475-600 mm**
- Ocean Ergo, Ergo XL and 24" Ergo model **Ocean 24" Ergo XL and VIP models
- No tools required

Accessories

Seating and backrests



Redesigned to mirror ergonomic seat plate, the softer feel gives users with sensory or skin protection issues added relief. Smaller aperture also available for shorter or thinner users. Easy to fit and clean.



The solid back support

Solid backrest

Soft seat insert

Available as a standard

Can be used with the soft

seat when the commode

function is not required.

Universal soft seat Comfortable, closed

seat surface for when offers an alternative to commode function is the standard tension not in use. adjustable backrest, when firmer support is required to achieve a more upright posture. Easy to clean.



Wider XL armrests Soft backrest cushion Increases the distance

For increased comfort. Can also be used to between the armrests by 40 mm on each reduce the overall seat width and depth, helping to position smaller users.



creates greater depth

of immersion promoting

increased pelvic stability

and positioning.

Armrest pads

New design offers

Comfort backrest

An alternative to the

a softer feel material

for increased comfort.

standard backrest with

improved support and

comfort for the user

(not for XL armrests).



XL back

Increases the width between the armrests by 80 mm. Can be retrofitted to standard Ocean Ergo models



Variable soft seat

Hygiene cut out can be positioned in all four directions, improving access for carers assisting with personal hygiene.



Safety



For additional safety and reassurance.



Safety lock for armrest Locks armrest into the vertical or horizontal



Front and rear anti-tippers For additional safety plus front and rear

stability.



Tilt-stop

An aid for those who are hoisted in/out of the chair, the tilt-stop screws ensure the frame remains horizontal.

Positioning



safety. Easy to attach. Machine washable.



Broad pad manufactured from soft material with single pull buckle. Can be or in addition across the chest for positioning and



Lateral supports

Adjustable in height via holder three possible positions on the backrest cane. Broad lateral pads can be repositioned for depth adjustment. Pad not included. comfort and effective pressure redistribution.



Ergonomic headrest Ergonomic headrest/ neckrest

Offers multidirectional Offers exceptional adjustment to correct management of the or accommodate for head position. more complex head/ Neckrest shown. neck positioning needs.

Toileting



with foot support Offers both height and

angle adjustment with a contoured pad for added positioning and comfort. Including a foot support for those who may need to elevate their lower limb(s) for medical reasons or comfort.



Calf and amputee support Calf and amputee



Splash guard Clips on easily for

use with the standard plastic seat, standard soft seat and small soft seat.



Toilet pan

Ergonomic handle allows the pan to be easily removed and replaced. Features an integrated lid holder, so the lid is no misplaced during use.



Ocean Ergo

Customise each of our **Ocean** Ergo chairs

Our Ocean Ergo range of shower chairs can all be customised to suit a wide number of needs. By making a simple adjustment or by adding one of our easy-to-fit accessories, each of our modular shower chairs are more versatile than ever before.



Ocean Ergo

An entry level attendant propelled model



Ocean Ergo XL

An entry level attendant propelled model with a heavier maximum user weight



Ocean 24" Ergo

An entry level self-propelling model



Ocean 24" Ergo XL

An entry level self-propelling model with a heavier maximum user weight



Ocean VIP Ergo

A tilt-in-space model for additional seating and positioning support



Ocean Dual VIP Ergo

A tilt-in-space model with recline for additional seating and positioning support

Ocean Ergo

Compatibility

- ✓ Compatible
- x Not compatible

Ocean Ocean

+ Included

Ocean

Ocean Ocean Ocean

Most accessori	es are cross compatible	è
with all Ocean	Ergo models.	

with all Ocean Ergo models.	Product code	Ergo	24" Ergo	Ergo XL	24" Ergo XL	VIP Ergo	Dual VIP Ergo
Standard soft seat	AP1603142	1	✓	✓	/	√	/
Soft seat insert	AP1603143	1	✓	√	✓	√	√
Small soft seat	AP1603145	1	✓	/	✓	√	/
Small soft seat insert	AP1603152	1	1	1	/	/	√
Splash guard	AP1603155	/	✓	✓	✓	√	✓
Universal soft seat	AP1603154	1	√	√	1	√	✓
Ergonomic special soft seat	AP60121311	/	✓	√	✓	√	✓
Variable soft seat	1544469	1	√	√	/	√	✓
Safety bar 2.0	AP60123732	/	✓	x	x	√	1
Safety bar 2.0 XL (only with armrests 2.0)	AP60123771	х	х	✓	√	Х	x
Safety bar XL (only for XL chairs with old armrests)	1470078	х	х	√	✓	Х	х
Chest belt	1470081	1	✓	✓	√	√	✓
Pelvic belt	1470082	1	✓	✓	1	√	√
XL back	AP1604015	1	√	+	+	х	x
Solid backrest	1575299	1	✓	X	х	√	√
Soft backrest cushion	1558769	1	✓	✓	✓	✓	✓
Armrest pad (1 pad; only for new armrests 2.0)	AP1643748	1	1	✓	1	✓	✓
Wider XL armrests	1470086	1	✓	✓	✓	✓	✓
Safety lock for XL armrest	1470077	1	✓	✓	✓	✓	✓
Safety lock 2.0 (for standard armrests 2.0)	AP1643760	1	✓	✓	✓	✓	✓
24" wheels (set)	16379	1	+	✓	+	X	x
Calf and amputee support right	AP1643573	1	✓	✓	✓	✓	✓
Calf and amputee support left	AP1643574	1	✓	✓	✓	✓	✓
Calf support with foot support right	AP1643543	1	✓	✓	✓	✓	✓
Calf support with foot support left	AP1643572	1	✓	✓	✓	✓	✓
Lateral support	1535077	1	✓	✓	✓	✓	✓
Ergonomic headrest holder	AP1531613	✓	✓	✓	✓	√	✓
Ergonomic headrest	1531465	1	✓	✓	✓	✓	✓
Ergonomic neckrest	1531608	1	✓	✓	✓	✓	✓
Toilet pan	10230	✓	✓	✓	✓	✓	✓
Toilet pan holder	AP1602078	1	✓	✓	✓	✓	✓
Front and back anti-tipper Ergo VIP models	AP1603992	x	x	x	x	✓	✓
Front anti-tipper Ocean Ergo models	AP1601986	1	✓	✓	✓	X	x
Comfort soft backrest VIP Ergo models	AP1603606	Х	Х	Х	х	√	√
Comfort soft backrest Ocean Ergo models	AP1604136	√	√	Х	х	X	Х
Comfort soft backrest Ocean Ergo XL models	AP1604016	Х	Х	✓	1	X	х
Tilt-stop screw set	AP60129156	х	х	х	х	1	✓

Tel: +44 61 487 70 70 hqeurope@invacare.com

hqeurope@invacare.com www.invacare.eu.com Invacare International GmbH Neuhofweg 51 4147 Aesch BL Switzerland

