

AQUATEC®

Ocean

Ergo Family



Modular personal care solutions

SAFE | INDEPENDENT | ERGONOMIC



Customer insight-led design

to enhance independence, safety and flexibility.

Personal care is highly private activity and is an important part of a person's life and wellbeing. The **Aquatec Ocean Ergo Family** of shower chair commodes have been designed with user comfort, independence and dignity in mind. Based on the well-established Aquatec Ocean range, the Ocean Ergo family comes packed with a range of new features, designed to make personal care activities easier for both the user and carer.

Safety, independence and dignity

A newly contoured seat surface has been developed to encourage a more upright and stable seating position, promoting user independence and overall

comfort. Enhanced stability at the front of the base allows for a smaller footprint and offers users added safety and reassurance when seated in the chair. Furthermore, the optimised wheelbase makes it even easier for a carer to manoeuvre the chair in small spaces as well as allowing for a better fit over standard toilets for independent toileting.

Modular and flexible design

The modular design of the Ocean Ergo family means that the chairs can be adapted to suit a wide range of user needs by making some simple adjustments, or by adding some of our easy to fit, off-the-shelf accessories.

OTHER features include:

- ▶ Stable, stainless steel frame
- ▶ Machine washable, tension adjustable back with quick-release buckles
- ▶ Easy disassembly for transportation and storage



NEW

FEATURES

1 Innovative ergonomic seat

- ▶ Insight-led design encourages a more upright seating position for enhanced independence and comfort

2 Greater adaptability

- ▶ With an increased height adjustment range on the footrests and armrests

4 Compact and easy to handle

- ▶ Optimised wheel base makes it easy to manoeuvre in small spaces

3 Extremely stable

- ▶ Enhanced stability at the front for additional safety and reassurance



NEW Contoured, ergonomic seat plate

- ▶ Based on the imprints people leave behind when seated
- ▶ A 5° seat tilt lifts the knees into a natural squat position for toileting
- ▶ Transfer supports at the front to allow the user to position themselves into a comfortable position



Highly Configurable

Available with a extensive range of positioning, safety, toileting and seating accessories that can be used to support a wide range of conditions and postural needs.

See full range of accessories online at www.invacare.eu.com



Positioning



Safety



Toileting



Seating and Backrests

Tilt in space models:

Ocean Ergo VIP and Ocean Ergo Dual VIP

- ▶ Both have a -5-40° seat tilt angle with a smart weight shift mechanism that ensures optimal stability on a small footprint, whilst the Ocean Dual VIP has an additional backrest recline of 5°-35°



Introducing the **Ocean Ergo** Family

From a standard transit model, through to tilt and recline, meet the flexible range of shower chairs that offer personal care solutions for a wide range of physical conditions.



Ocean 24” Ergo XL ▲

An entry level, plus-size model with a maximum user weight of 180 kg

Ocean Ergo XL ▲

An entry level plus-size attendant propelled model, with a maximum user weight of 180 kg

Ocean VIP Ergo ▲

Designed for individuals who require slightly more assistance with postural control when seated

Ocean Ergo ▲

An entry level attendant propelled model

Ocean Dual VIP Ergo ▲

Designed for individuals with complex seating and positioning needs

Ocean 24” Ergo ▲

An entry level self-propelling model

Technical data ▼

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.



Seat width ▼



Seat depth ▼



Total width ▼



Total depth ▼



Width between armrests ▼

Ocean Ergo	480 mm	480 mm	565 mm	905 mm	455 mm
Ocean Ergo XL	480 mm	480 mm	645 mm	905 mm	535 mm
Ocean 24" Ergo	480 mm	480 mm	670 mm	1050 mm	455 mm
Ocean 24" Ergo XL	480 mm	480 mm	670 mm	1050 mm	535 mm
Ocean VIP Ergo	480 mm	480 mm	565 mm	1000 mm	455 mm
Ocean Dual VIP Ergo	480 mm	480 mm	565 mm	1000 mm	455 mm



Tilt angle seat ▼



Tilt angle back ▼



Weight ▼



Load capacity ▼

Ocean Ergo	-	-	17 kg	150 kg
Ocean Ergo XL	-	-	18 kg	180 kg
Ocean 24" Ergo	-	-	23 kg	150 kg
Ocean 24" Ergo XL	-	-	24 kg	180 kg
Ocean VIP Ergo	-5° - 40°	-	23.5 kg	150 kg
Ocean Dual VIP Ergo	-5° - 40°	5° - 35°	28.3 kg	150 kg

Frame colour ▼



White

Regulatory data ▼



Tel: +41 61 487 70 70
hqeuropa@invacare.com
www.invacare.eu.com

Invacare International GmbH
 Benkenstrasse 260
 4108 Witterswil
 Switzerland

