

Propad™ Low Profile Cushion



Profiled foam cushion designed
for comfort and support



Yes, you can.®

The *Invacare Propad Low Profile cushion* provides users with an option of having a lower cushion, rather than a standard height cushion. Some users prefer a low profile cushion, finding it more compatible with their daily living activities.

Features and Options



Cover helps to prevent a build up of shear and friction forces against the user

Cover is two-way stretch, water-resistant and vapour permeable. Cover features an anti-slip base.

Available in a range of sizes

The *Invacare Propad Low Profile Cushion* is available in 102 sizes, ranging from 250 x 250 mm – 600 x 550 mm.

Full length flap conceals the zip

The full length flap concealing the zip reduces the potential for fluid ingress into the cushion and protects the foam core. The cover can be removed and laundered.

Castellated Foam Cushion




Each castellation delivers enhanced support and pressure area care to the user.

Reduces shear and friction forces

The castellations in the *Invacare Propad Low Profile cushion* reduce shear and friction forces on tissue.

Technical data


For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

<p>Size range</p>  <p>A range of 102 sizes</p>	<p>Max. user weight</p>  <p>76 kg</p>	<p>Max. cleaning temperature</p>  <p>80°</p>
---	--	---

Propad Low Profile Cushion

Manufactured to European Standards:
BS5852: Crib 5
BS7175: Crib 5

Colours



Tel: +44 61 487 70 70
hqeurope@invacare.com
www.invacare.eu.com

Invacare International GmbH
Neuhofweg 51
4147 Aesch BL
Switzerland



Yes, you can.®