

Propad™ Original Cushion



Profiled foam cushion designed
for comfort and support



The *Invacare Propad Original Cushion* has been designed to facilitate good pressure redistribution when used with both wheelchairs and standard armchairs. This cushion is constructed to ensure improved support and prolonged longevity.



Features and Options



Cover helps prevent a build up of shear and friction forces against the user

Cover is two-way stretch, water-resistant and vapour-permeable. Cover features an anti-slip base.



Available in a range of sizes

The *Invacare Propad Original Cushion* is available in 3 standard sizes, 410 x 410, 430 x 430, and 460 x 410 mm. (See size chart for details).



Full length flap conceals zip

The full length flap concealing the zip reduces the potential for fluid ingress into the cushion and protects the foam core. The cover can be removed and laundered.

Deep contoured profiled foam cushion designed to deliver enhanced pressure care

The *Invacare Propad Original Cushion* is made from a layer of castellated foam, each castellation delivering enhanced support and pressure area care to the user. The “keyhole” cut at the base of each castellation ensures good air circulation, reducing potential for skin maceration.

Contoured seating surface

Offers improved stability and postural support; castellations reduce shear and friction forces on tissue.

Technical data

For more comprehensive pre-sales information about this product, including the product’s user manual, please see your local Invacare website.

	Size range	Height	Total product weight	Max. user weight	Max. cleaning temperature
					
Propad Original Cushion	406 x 406 mm 432 x 432 mm 406 x 457 mm	102 mm	1 kg	108 kg	80°

Manufactured to European Standards:
Interior Foam: BS EN 597-1 Cig, BS EN597-2 Match, BS 5852 (Crib 5)
Cover: BS 7175 (Crib 5)

Colours

Blue PU

Black Towelling



QR Code

Tel: +41 61 487 70 70
hqeurope@invacare.com
www.invacare.eu.com

Invacare International GmbH
Benkenstrasse 260
4108 Witterswil
Switzerland



Yes, you can.®