

# **Propad™ Profile Cushion**



Profiled foam cushion designed  
for comfort and support



Yes, you can.

The *Invacare Propad Profile cushion* provides users with an option of having a lower 3" high cushion, rather than a standard height (4") cushion. Some users prefer a lower profile cushion, finding it more compatible with their daily living activities.


### Castellated Foam Cushion

Each castellation delivers enhanced support and pressure area care to the user.

### Reduces shear and friction forces

The castellations in the *Invacare Propad Profile cushion* reduce shear and friction forces on tissue.

## Features and Options



### Cover helps to prevent a build up of shear and friction forces against the user

Cover is two-way stretch, water-resistant and vapour permeable. Cover features an anti-slip base.



### Available in a range of sizes

The *Invacare Propad Profile Cushion* is available in 102 sizes, ranging from 25 x 25 cm (250 x 250 mm) – 60 x 55 cm (600 x 550 mm).






### Full length flap conceals the zip

The full length flap concealing the zip reduces the potential for fluid ingress into the cushion and protects the foam core. The cover can be removed and laundered.


## Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local *Invacare* website.

	Size range	Max. user weight	Max. cleaning temperature
<b>Propad Profile Cushion</b>	 A range of 102 sizes	 95 kg	 80°

Manufactured to European Standards: BS5852: Crib 5  
BS7175: Crib 5

### Colours



Tel: +44 61 487 70 70  
hqeurope@invacare.com  
www.invacare.eu.com

Invacare International GmbH  
Neuhofweg 51  
4147 Aesch BL  
Switzerland



Yes, you can.®