



Invacare®

XLT Max

Stability, Rigidity, Security: Be Active!



The new *Invacare XLT Max* is a bariatric, rigid active wheelchair ideal for day-to-day use. Built with a titanium frame, it is designed to support a user weight of max. 180 kg with seat width possibilities of up to 60 cm. The *XLT Max* offers users various levels of support thus enabling greater mobility and a more active lifestyle.

The Perfect Fit

Measured and constructed to fit the requirements of each user, the *Invacare XLT Max* is an ideal chair for those previously unable to take advantage of an active wheelchair. Built for maximum security and support in combination with good driving performance and manoeuvrability, the *XLT Max* is the only active rigid bariatric wheelchair on the market today.



Simply Compatible

The *XLT Max* is a rigid titanium wheelchair offering users stability and the utmost rigidity. The well-designed fit of the wheelchair helps users with a weight of up to 180 kg retain and improve mobility and independence in their daily lives. Users will benefit from the extensive and readily available array of Invacare options which can easily be fitted to the *XLT Max*.

Features and Options



Invacare Flo-tech Contour

Ideal bariatric cushion solution. Improved posture and increased support surface area significantly contributes to preventing pressure ulcers.



Footplate converter

Divided into three parts, the footplate comes with a removable mid section, to facilitate transfers.

Anti-tipper device

The anti-tipper device is height adjustable and can also act as a step tube.




Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

Invacare XLT Max	Seat width	Seat depth	Seat angle	Seat height	Backrest angle	Backrest height	Length of footrests	Driving unit width
	505 / 555 / 605 mm	500 / 550 / 600 mm	0 - 14°	440 - 520 mm	± 10°	400 - 490 mm	380 - 480 mm 270 - 370 mm high mounted	Seat width + 210 mm
Invacare XLT Max	Total length inc. Footrests	Total height	Max. User weight	Total weight	Crash test			
	930 - 1080 mm	790 - 1100 mm	180 kg	15,5 kg	Crash-tested ISO 7176-19			

Frame colour



QR Code

Tel: +41 61 487 70 70
hqeurope@invacare.com
www.invacare.eu.com

Invacare International GmbH
Benkenstrasse 260
4108 Witterswil
Switzerland

