

# True Track<sup>®</sup> Plus



... for those who want more!



The new generation **True Track Plus** motor gives more driving range and power for drivers who enjoy a challenge. **True Track Plus** immediately responds to the user's commands and provides more freedom to discover the world by easily and smoothly overcoming obstacles. Users will appreciate the speed and discretion that comes with this absolutely silent motor.

### Unrivalled efficiency, power and torque

Improved efficiency gives more range and power. Easier slope climbing at slow speeds with up to 40% more torque. Power and immediate response are maximised for beginners or users of alternative driver controls as well as for the more experienced driver who enjoys a challenging, dynamic and fast ride.

### Absolutely silent and maintenance free

No friction, no wearing parts, no carbon brushes – no maintenance. This option completely removes the noise and vibrations that occur on conventional geared motors. The gearless, brushless motor provides maximum discretion for the customer.



## Features and options



### Speed choice

Adapt to your client's evolution and easily programme all speed choices without loss of torque.



### Smarter design

Compact motor is hidden by the drive wheel it attaches to.



### ACS2 controls

**ACS2** controls offers a number of input options.

**True Track Plus** is available for *Invacare Storm®4 X-plore* and *Invacare Storm®4*.

## Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local **Invacare website**.

Motor capacity	Speed	Driving range

**True Track Plus**      600 W      6 / 10 / 13 km/h      up to 59 km<sup>1</sup>

1. According to ISO 7176-4 on **Storm4** and **Storm4 X-plore** programmed to 6 km/h



Scan here to access the website.

Tel: +41 61 487 70 70  
[hqeuropa@invacare.com](mailto:hqeuropa@invacare.com)  
[www.invacare.eu.com](http://www.invacare.eu.com)

Invacare International GmbH  
 Benkenstrasse 260  
 4108 Witterswil  
 Switzerland

