

Dolomite[®] Brass



Lightweight and reliable for everyday use

COMPACT | EFFORTLESS MANOEUVRABILITY | ADJUSTABLE



Yes, you can.™

Dolomite Brass

A lightweight aluminium rollator that ensures effortless manoeuvrability over various surfaces both indoors and out, the **Dolomite Brass** is the perfect companion for every day use. Designed with the user in mind, the soft padded seat provides a sturdy surface for those who need a rest whilst it can be easily flipped up when folded for easy storage.

Compact for those who are always on the move

Able to stand up-right and weighing only 7kg*, the '**Brass**' can be easily lifted for transport or conveniently stored in compact spaces.

Robustly designed for peace of mind, the strong aluminium frame has been extensively tested for every day use and is also fitted with puncture proof tyres.

Suitable for a variety of individuals

An impressive maximum user weight of 135kg and height adjustable handles ensure that a variety of user sizes can be comfortably accommodated.



Standard features



Crutch holder

Makes the crutch easy to access without getting in the way



Soft fabric basket

Enables users to carry small personal items during use



Kerb climber

Helps to simplify passing kerbs or similar elevated obstacles

KEY **EXCLUSIVE**
FEATURES




Technical data ▼

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

	Width between outside wheels ▼	Depth ▼	Width between wheels ▼	Overall height ▼	Seat height ▼	Width between handles ▼
Dolomite Brass	640mm	620mm	580mm	815 - 995mm	566mm	486mm
	Wheel diameter ▼	Product weight* ▼	Max. user weight ▼	Turning diameter ▼	Diameter of frame tube ▼	Wheel dimensions ▼
Dolomite Brass	191mm	7kgs	135 kg	817mm	Ø33×18 x 1.2	200mm

*Weight excludes any fitted accessories

Frame colours ▼



RAL 5022

Please note: colours may vary slightly from those displayed here.

Tel: +41 61 487 70 70
hqeuropa@invacare.com
www.invacare.eu.com

Invacare International GmbH
 Benkenstrasse 260
 4108 Witterswil
 Switzerland

© 2019 Invacare International GmbH. All rights reserved. All information quoted is believed to be correct at time of print.
 Invacare reserves the right to alter product specifications without prior consultation. Brass - EU - 01/2020



Yes, you can.®

