

Perfect combination of style and strength



LIGHTWEIGHT | VERSATILE | MODERN

# Walk the walk with Gloss

**Dolomite** is advancing its range of established rollators by introducing the lightweight **Gloss** family. Designed to support users while out and about, the **Dolomite Gloss** is the ideal rollator to empower independence. Stable and robust, the rollators contain all the core values of the established **Dolomite** range whilst encompassing a modern and sleek appearance.

Available in 4 seat heights ranging from 450mm – 680mm, the **Dolomite Gloss** can accommodate a variety of clients up to 160kg. An extensive range of optional accessories are available to quickly adapt the versatile rollator to fit a range of conditions and needs.

Complementing the 160kg maximum user weight is an extremely lightweight and agile 6.5kg aluminium frame\* frame that makes transportation and storage effortless. The attractive cross-folding frame offers a generous walking space, while the firm patented seat provides comfort for clients that require extra support. For additional safety, the **Dolomite Gloss** has built-in audible and visual feedback to accompany the folding and unfolding mechanism.

With its versatility, quality and stability, the **Dolomite Gloss** will empower all users to walk the walk in safety and style.







#### Easy transportation and storage



The frame of the *Gloss* folds to an extremely narrow 190mm<sup>\*\*</sup>, making it perfect for storing in compact spaces. The design of the rollator means that it can also be stored securely in an Upright position when not in use.

For further independence, the lightweight of the *Gloss* ensures that it can be easily lifted into the boot of a car by client or carer for transportation between locations.

\*Gloss 600 weights 6.5kg. \*\*520, 600 & 680 size models



#### Height adjustable, ergonomic handles

The forward facing, ergonomically designed handles can be height and angle adjusted to accommodate for a variety of client sizes.

# Large walking

space

needs.

 Comfortable walking space inside the frame promotes strong posture and freedom during movement.



 Patented 3-piece seat design facilitates simple and compact folding and unfolding

# Lightweight cross folding design

Weighing only 6.5kg and measuring only 190mm when folded\*, the rollator can be easily lifted and stored in an upright position.

\*600 size model

### Additional **KEY** features

Versatility for individual needs

Available in 4 different seat heights, while an

extensive range of accessories ensures that the

rollator can be adapted to a client's condition or

- Unobtrusive wiring enhances sleek and modern appearance
- The combination of dynamic and parking brakes makes stopping and resting effortless
- Safety locking strap for secure storage when not in use
- Easy to reach foldable storage basket available for transportation of goods
- Enhanced safety and independence due to numerous accessories.
- Cross compatibility of components with future models within the *Gloss* family





# 4 key features and benefits

#### Lightweight cross folding design

Alongside being extremely light, the **Gloss** has an impressive maximum user weight of 160kg. It's small cross folding design ensures the frame is very narrow when folded upright, enabling the rollator to be easily stored in tight spaces.

The lightweight design also enables clients and carers to easily lift the *Gloss* into the boot of a car for transportation.



#### Versatility for individual needs

The *Gloss* is available in 4 seat heights: 450mm, 520mm, 600mm and 680mm. These varying seat heights ensure that the *Gloss* can be suitable for a range of client sizes. The versatility of the rollator is further enhanced by an extensive range of accessories that can quickly adapt the frame depending on a client's condition or specific needs.



#### Large walking space

A signum for Dolomite rollators, the **Gloss** has a large walking space inside the footprint of the frame. Designed to promote strong posture, the comfortable area helps to provide users with support and freedom whilst walking.



#### Height adjustable, ergonomic handles

The height of the handles on the *Gloss* can quickly and easily be adjusted to satisfy a range of heights by merely unscrewing the knob positioned on the outside of the frame.

The ergonomic, forward facing handles can also rotate inwards 20° to accommodate a range of user needs.



# Introducing the **Dolomite Gloss** family

#### The new rollator in *Dolomite*'s esteemed range. Designed for style and strength.



#### Available options and accessories:

- Back strap: Comfortable, safe and supportive whilst sitting.
- Cane holder: For ease of use for those with additional needs.
- Foldable storage basket: Can easily be removed or be kept on the rollator while folded.
- Tray: Detachable tray for transportation of food and drinks to other surfaces.
- Kerb climber: Additional support while climbing kerbs or uneven surfaces.
- Hemiplegic handle set: Optional adaptation for one-handed or stroke clients.

- Slow down brake: Provides extra support for when braking is challenging.
- Brake booster: Provides extra braking force for those with weak hands.
- Soft padding: Additional seat padding for extra comfort.
- Oxygen Tank Holder: Allows for secure mounting of oxygen tanks.
- Crossbar: Makes it possible to handle the rollator with one hand.
- One Hand Brake: Makes it possible to activate both brakes with one hand.
- IV Pole: Securely fixes IV drip.





 Hard basket: Durable basket available for the 520, 600 & 680 models.

## Technical data -

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

E		E
Width between outside wheels ▼	Depth 🔻	Width b whee







	Width between outside wheels 🔻	Depth 🔻	Width between wheels ▼	Overall height 🔻	Seat height 🔻
Dolomite Gloss 450	518 mm	615 mm	423 mm	635 - 770 mm	450 mm
Dolomite Gloss 520	595 mm	675 mm	501 mm	725 - 860 mm	520 mm
Dolomite Gloss 600	602 mm	720 mm	509 mm	795 - 920 mm	600 mm
Dolomite Gloss 680	609 mm	790 mm	501 mm	915 - 1050 mm	680 mm









	Seat width 🔻	Wheel diameter 🔻	Max. user weight 🔻	Product weight* ▼
Dolomite Gloss 450	370 mm	200 mm	160 kg	5.9 kg
Dolomite Gloss 520	445 mm	200 mm	160 kg	6.3 kg
Dolomite Gloss 600	445 mm	200 mm	160 kg	6.5 kg
Dolomite Gloss 680	445 mm	200 mm	160 kg	6.8 kg

\*Excludes any accessories

# Frame colour ,



Tel: +41 61 487 70 70 hqeurope@invacare.com www.invacare.eu.com

Invacare International GmbH Benkenstrasse 260 4108 Witterswil Switzerland



© 2019 Invacare International GmbH. All rights reserved. All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation. Dolomite Gloss - EU - 08/2019