

AQUATEC® PICO GREEN



100% Natural antibacterial effect

ANTIBACTERIAL | COMFORTABLE | SUSTAINABLE



The *Aquatec* **PICO GREEN** shower chair is a pioneering product in shower aids, available as standard version or as a commode including toilet pan. With a modern and unique design, it can complement the most elegant of bathrooms whilst providing excellent comfort for those who need a little extra safety and support whilst showering. Lightweight yet strong enough to hold up to 160 kg, the **PICO GREEN** can be adjusted to suit a range of heights and seating positions.



Kind to the environment, kind to your skin

The real advantage of the **PICO GREEN** is that its material is made from 10% natural wood fibres sourced from sustainable forests in Germany. Not only does this provide an initial fragrant aroma but offers a 100% natural antibacterial effect to minimise the risk of infection and reduce the need for strong chemical disinfectants.

The wooden fibres have a protective function and kills bacteria without utilising any chemical additives or nano-silver ions.

The natural antibacterial properties of the **PICO GREEN** have been proven by an independent biology laboratory to kill 99.9% of bacteria*!

*measured on E. coli and S. aureus



Technical data ▼

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.



Seat height ▼



Seat width ▼



Seat depth ▼



Total width ▼



Total depth ▼



Width between armrests ▼

**PICO GREEN /
PICO 3in1 GREEN**

425 - 575 mm

430 mm

420 mm

575 mm

520 mm

430 mm



Weight shower chair ▼



Max. user weight ▼



Colours ▼



White with brown wooden fibres

PICO GREEN

4.6 kg

PICO 3in1 GREEN

5.4 kg

160 kg

Tel: +41 61 487 70 70
hqeurope@invacare.com
www.invacare.eu.com

Invacare International GmbH
Benkenstrasse 260
4108 Witterswil
Switzerland



Yes, you can.®