

**rea**<sup>®</sup>

# Azalea<sup>®</sup> Tall



No limits on comfort for  
taller users



Yes, you can.

**Rea Azalea Tall** is specially designed to meet the needs of tall users who require a Tilt in Space wheelchair with a longer seat support. Adapted from the **Rea Azalea**, the **Rea Azalea Tall** boasts all the advantages of a reliable, tilting wheelchair and offers a unique weight-shifting mechanism. The **Rea Azalea Tall**: at the top of its class for comfort and stability. The improved Flex 3 backrest adds individual adaptability. 5 frame colours and three upholsteries are available.

## Real compatibility

The **Rea Azalea Tall** has all the features of the **Rea Azalea** such as Tilt in Space weight-shifting for a really stable tilt. All accessories such as rear wheels for self-propelling or smaller wheels for transport and a different seat system also fit the **Rea Azalea Tall**.

## Easy maneuverability

The adjustable push bar has been designed with the carer in mind, it offers a good ergonomic grip and can be adjusted to suit the needs of a carer.



### Longer seat plate

An extra 50 mm of seat depth is designed to support longer legs.



### Chassis length

An extra 50 mm in length enables greater forward stability.



### Seat height

An extra 50 mm of height allows for longer leg-rests.

For other Rea Azalea accessories, please consult the Rea Azalea brochure or visit [www.invacarearea.com](http://www.invacarearea.com) to see videos of them

## Technical data


For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

|                    | Seat width ▼   | Seat depth ▼   | Seat height ▼                 | Backrest height ▼  | Armrest height ▼   | Length of footrests ▼ | Seat angle ▼               | Backrest angle ▼                         |
|--------------------|--|----------------|-------------------------------|--|--------------------|-----------------------|----------------------------|--|
| <b>Azalea Tall</b> | 390-450/ 440-500/<br>490-550 mm<br>(390-490/ 440-540/<br>490-590 mm) | 480 - 550 mm   | 500 mm                        | Backrest plate:<br>560 - 790 mm<br>Tension adjustable<br>backrest:<br>540 - 715 mm | 240 - 360 mm       | 330 - 627 mm          | -1° to 23°<br>(+1° to 25°) | 0° - 30°                                 |
|                    |  |                |                               |  |                    |                       |                            |  |
|                    | Driving unit width ▼   | Total height ▼ | Total length inc. Footrests ▼ | Total weight ▼   | Max. User weight ▼ | Transport weight ▼    | Crash tested ▼             | Maximum slope for use of parking brake ▼ |
| <b>Azalea Tall</b> | SW + 220 mm<br>(SW + 230 mm)   | 1000 - 1300 mm | 950 - 1085 mm                 | 34 kg<br>(SW 440 mm)   | Max 135 kg         | 18,7 kg <sup>1</sup>  | ISO 7176-19                | 7°                                       |

## Frame colour



## Upholstery



1. Transport weight increases with 4,2 Kg when Power tilt & recline option is used



Tel: +44 61 487 70 70  
[hqeuropa@invacare.com](mailto:hqeuropa@invacare.com)  
[www.invacare.eu.com](http://www.invacare.eu.com)

Invacare International GmbH  
 Neuhofweg 51  
 4147 Aesch BL  
 Switzerland

Rea, Rea design and DSS (Dual Stability System) design are registered trademarks of Invacare International. The **Rea Azalea** is designed with DSS (Dual Stability system). It provides the user with enhanced stability and safety when being tilted or the backrest is reclined.

