

# Toilet seat raisers

a range of adaptable toileting solutions

AQUATEC

#### Providing dependable toileting support



Adjustable, reliable, comfortable

# **Introducing the Aquatec range** of raised toilet seats

#### An adaptable range offering supportive toileting solutions.

Providing a reliable and supportive toileting solution for individuals is very important to aid independence.

The dependable Aquatec range offers a secure and stable environment with fixations that support the placement of the raised toilet seats in situ.

The range offers a variety of possible seat heights to accommodate individual sizes and to ensure optimal comfort.


Featuring models that offer supportive armrests for stability and standing transfers, including the swing-away options for ease of functionality.

### Aquatec 90000

AQUATEC Toilet seat raisers

- > Available in three different seat heights: 2 cm, 6 cm, 10 cm.
- > Can be ordered with or without armrests.
- > Swing-away, ergonomic armrests for additional support when sitting down and standing up.
- > Innovative clip system: The lid and toilet seat can be removed from the base; no tools are required.







# **Technical data**

For more comprehensive pre-sales information about this product range, including the product's user manuals, please see your local Invacare website.

	$\mathbf{H}$	F		
	Seat width 🔻	Seat depth 🔻	Seat height 🔻	Total width 🔻
Aquatec 900	360 mm	425 mm	60 / 100 / 150 mm	590 mm (with armrests) 430 mm (without armrests)
Aquatec <b>90000</b>	365 mm	390 mm	20 / 60 / 100 mm	540 mm (with armrests) 380 mm (without armrests)

				$\bigcirc$
	Total depth 🔻	Total height 🔻	Width between armrests 🔻	Width of opening 🔻
Aquatec 900	490 mm	380 mm (with armrests) 180 mm (without armrests) 200 mm (seat to armrests)	510 mm	225 mm
Aquatec <b>90000</b>	475 mm	230 / 270 / 310 mm (with armrests) 30 / 70 / 110 mm (without armrests) 200 mm (seat to armrests)	450 mm	205 mm

	Depth of opening 🔻	Weight 🔻	Load capacity 🔻
Aquatec 900	285 mm	5.0 kg (with armrests) 2.7 kg (without armrests)	120 kg
Aquatec <b>90000</b>	275 mm	4.4 kg (SH:100 mm) 4 kg (SH:60 mm) 3.6 kg (SH:20 mm) (with armrests)	150 kg

hite		
	hite	hite

Frame colour

#### Tel: +41 61 487 70 70 hqeurope@invacare.com www.invacare.eu.com

Invacare International GmbH Neuhofweg 51 4147 Aesch BL Switzerland

 $\odot$  2024 Invacare International GmbH. All rights reserved. All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation. EU – 04/24



Yes, you can: