

Standard

Option field

**Required fields:**<sup>1</sup>

Information required for processing

ORDER

QUOTE

**Invoice address** (stamp or fill in legibly)

Cust. No.

Company

Street

ZIP/Place

**Shipping address** (if different from invoice address)

Company

Street

ZIP/Place

### Type of order

Self-fitting<sup>2</sup> (fully complete the section "Information about the wheelchair"!)

Mounting at the Alber factory<sup>3</sup> 1592525 **100,00 EUR**

**e-motion DuoDrive**  
With control unit (bracket included) 1594352 **please refer to your price list**  
Wheel size 22" (550 mm)

**e-motion DuoDrive**  
With control unit (bracket included) 1594353 **please refer to your price list**  
Wheel size 24" (600 mm)

**e-motion DuoDrive**  
With control unit (bracket included) 1594355 **please refer to your price list**  
Wheel size 25"  
(625 with Right Run tyres, 615 with Marathon Plus tyres)

**Upcharge ECS Remote Control** 1592486 **500,00 EUR**  
Rollback delay can only be activated with  
ECS Remote Control or smartphone app with Mobility Plus Package

### Important accessories

Bracket for mounting to the wheelchair 1467178 **205,00 EUR**

Anti-tippers with jack-up function 1593206 (set of two) **265,00 EUR**

### Wheelchair details

Wheelchair manufacturer (e.g. Meyra)

Model (e.g. Eurochair, 1.850)

Year of manufacture

Version (e.g. Hemi)

Drum brakes<sup>5</sup>  Yes  No  
↳ Attention: Drum brakes must be removed

Car restraint system available (please consult with Alber customer service)

Wheel base extension available

Existing wheel position

1 = back

2 = centre

3 = front



**Delivery**  Standard 0,00 EUR

**Consignment**

Extended guarantee Li-Ion battery pack 1592523 **95,00 EUR**  
(to 36 months)

Seat width (seat tube measured externally)   cm

Seat height back (measured to lower edge of seat sling)   cm

Diameter footrest for bracket of control unit DuoDrive   mm

Wheel camber (if possible)

0°

3°

**Comments** (mounting, delivery etc.)

Standard

Option field

**Required fields:**<sup>1</sup>

Information required for processing

### Tyres

- Schwalbe Right Run (Standard) black 1592530 0,00 EUR  
(1" for wheel size 22" and 24", 1,1" for wheel size 25")
- Schwalbe Marathon Plus Evolution 1565894 48,00 EUR  
with SmartGuard puncture protection, black  
(not available for 22")
- Pro Speed puncture-proof solid PU 1" 1592532 125,00 EUR  
made of solid PU, black  
(+ 300 g per wheel, not available for 22" and 25")
- Puncture-proof with inlay 1 3/8" 1592533 125,00 EUR  
with PU insert and tyre cover, schwarz  
(+ 500 g per wheel, not available for 25")

### Push rims

- Stainless steel (standard) 1592526 0,00 EUR
- Coated push rim Tetra 1592527 175,00 EUR  
(not available for 25")  
(+250 g per wheel)
- Carbolife CURVE L 1592528 380,00 EUR  
powder coated, black  
(not available for 22")
- Carbolife QUADRO 1592529 460,00 EUR  
Tetra-Gripp coating  
(Only available for 24")

### Push rim covers (retrofit)

- push rim cover 22" 1592425 (pair) 125,00 EUR
- push rim cover 24" 1592421 (pair) 125,00 EUR
- push rim cover 25" 1592426 (pair) 125,00 EUR

### Spoke cover (retrofit)

- Spoke cover 22", Design „Comic“ 1592266 (pair) 145,00 EUR
- Spoke cover 24", Design „Sphere“ 1592268 (pair) 145,00 EUR
- Spoke cover 25", Design „Sphere“ 1592269 (pair) 145,00 EUR
- Spoke cover 24", plain 1592553 (pair) 145,00 EUR

### Accessories

- Ergonomic stub axle 1592972 (pair) 105,00 EUR
- Transport Axle 1591362 50,00 EUR
- Conversion of manual wheels to 1/2" stub axle<sup>4</sup> 1565005 85,00 EUR
- Brake extension kit 1467780 (pair) 22,00 EUR

Date  /  /  Signature of ordering party

Date 01.05.2021. All prices before VAT, including packaging ex factory. The general Terms and Conditions of Alber GmbH apply. Subject to changes.

1 In order to smoothly and quickly process your order we require ALL information in the fields highlighted in blue.  
 2 Installation kits can only be shipped to trained specialist dealers (installation training).  
 3 Please send the wheelchair with a copy of this order as soon as possible and freight paid to Alber.  
 4 If the manual wheels are used in combination with e-fix, they must be equipped with 1/2" stub axles, without drum brakes and tyres of size 1".  
 5 If present, drum brakes must be dismantled as they are not compatible with e-motion drive wheels. The use of standard lever brakes is required.