



AQUATEC®

90 Ergo

a new ergonomic raised toilet seat

Designed for comfort and ease of use

A discreet, modern style with no fixtures or fittings



Yes, you can.®

AQUATEC®

90 Ergo

The newest addition to our Aquatec range

Our obsession with quality means all Aquatec designs focus on the details. Our range has gained an unsurpassed international reputation. Every Aquatec solution is carefully designed to overcome personal care challenges.

Discover the new Aquatec 90 Ergo.



A new raised toilet seat with unrivalled features

The Aquatec 90 Ergo is our newest addition to our Aquatec range. It features an optimal 10cm seat height, new ergonomic design and is styled for both modern and traditional toilet designs making it extremely versatile. It strikes the perfect balance between cost effectiveness and a quality design. A new generation of raised toilet seat has arrived.



A focus on ergonomics

- > Shaped for comfort
- > Contours around the body
- > 5° tilt for easier toileting



ERGONOMIC
DESIGN

5°



Styled without compromise

The days of choosing between a functional or a stylish bathroom are over. Welcome to a new era of toileting accessories designed to promote independence without compromising on bathroom aesthetics.



Ergonomically designed

Aquatec 90 Ergo features an anatomically designed seat also found within our popular Aquatec Ocean Ergo shower chair range. The ergonomic seat perfectly shapes to the body's contours ensuring maximum stability and higher levels of comfort. The 5° tilt also promotes an improved position for easier toileting.



Optional lid

The Aquatec 90 Ergo is available with an optional lid for added aesthetic appeal and hygiene purposes.



Optional lid



Fits a range of toilet designs

The Aquatec 90 Ergo fits both the more traditional round as well as today's modern square-shaped toilets, without the need for clips or brackets.



Eco-conscious

Making more eco-friendly choices means being kinder to our environment as well as to our home. The Aquatec 90 Ergo is made with 47% recycled material for a reduced CO₂ footprint. Choosing Aquatec 90 Ergo means conserving resources for a more sustainable environment.



No fixtures or fittings

Designed with new friction-gripping, the Aquatec 90 Ergo ensures a highly secure and stable position without brackets, fixtures and fittings. Simply lift off and store when not in use.



No fixtures
or fittings



Compact for
easy storage








A reputation for quality and innovative design

Aquatec has long been synonymous with premium quality equipment for assisted and independent personal care. The range includes toileting, showering and bathing solutions. You can explore the full range by visiting www.invacare.eu.com.

Technical data

For more comprehensive information about this product, including the user manual, please visit your local Invacare website.

90 Ergo

				<div>White</div> <div>Colour</div>
Width of opening ▼	Total width ▼	Depth of opening ▼	Total depth ▼	
230 mm	410 mm	315 mm	411 mm	
				
Raised seat height (from toilet bowl) ▼	Total product height ▼	Weight ▼	Load capacity ▼	
100 mm	143 mm	1.6kg (with lid) 1.2kg (without lid)	225kg	

Tel: +41 61 487 70 70
hqeurope@invacare.com
www.invacare.eu.com

Invacare International GmbH
Neuhofweg 51
4147 Aesch BL
Switzerland



Yes, you can.®