


ULTRA-LIGHTWEIGHT, EASY TO USE AND COMFORTABLE



- ▲ Air valve automatically adapts the cushion to your body shape and weight
- ▲ Two compartment layout limits sliding
- ▲ Easy to clean and disinfect

The Vicair Liberty is an easy to use and comfortable wheelchair cushion. It is designed for users who require comfort in combination with a low to medium level of support and pressure injury prevention. The cushion is ready-for-use, rain proof and with two compartments it provides an even pressure distribution. This makes it an ideal cushion in a large variety of situations.

The automatic valve makes sure that the cushion adjust to every seating plane and every body contour.

This Vicair cushion is extremely lightweight.

PRODUCT PERFORMANCE

Skin protection / Anti decubitus

Pressure redistribution	▲▲▲△△
Shear force reduction	▲▲△△△
Microclimate control	▲▲△△△
Impact & vibration damping	▲▲▲△△

Positioning

Lateral stability	▲▲△△△
Frontal stability	▲▲▲△△
Sliding reduction	▲▲▲△△

Allows manual adjustments N/A

Increasing endurance ▲▲▲△△

Maintenance & Safety

Durability	▲▲▲△△
Reliability	▲▲▲△△
Functionality when misused	▲▲▲△△
Easy to maintain	▲▲▲▲▲
Level of hygiene	▲▲▲▲△

Functionality & Application

Comfort	▲▲▲▲△
Ultra light weight	▲▲▲▲▲
Allows easy transfers	▲▲▲△△
Available standard sizes	▲▲▲△△
Maximum user weight	▲▲▲△△

TECHNICAL SPECIFICATIONS

Thickness	5 cm
Average weight	350 grams
Max. user weight	150 kg

CUSHION SIZES

- Standard
- Not available

Inch	width										
	14	15	16	17	18	19	20	21	22	23	24
CM	35	38	40	43	45	48	50	53	55	58	60
14	35										
15	38										
16	40										
17	43										
18	45										
19	48										
20	50										
21	53										
22	55										
23	58										
24	60										

AVAILABLE COVERS

Comfair Cover and Incotec Cover.