

Action[®] 2_{NG}

TRANSIT LITE



Lighter and more compact

Even better for travelling



Yes, you can.®

Take *Action*[®] and Go.

The new Action 2 NG Transit Lite is the newest chair from our award-winning Action range. Comfortable and reliable, Action chairs open a world of potential. Each chair in the range has a highly versatile core design with customisable options. The Action range is recognised simply as one you can count on.



Say hello to the new Action 2 NG Transit Lite

The new Action 2 NG Transit Lite is now lighter, more compact and super easy to handle. Its reduced frame size makes it even better for travelling and fitting into snug storage spaces.

Cross compatible with the range of Action 2 NG options, it's easy to adjust and maintain to suit your needs. Supporting you wherever you need to go.

Let's embrace new adventures.



1



Lighter and more compact.

We've reduced the weight and frame size, making Action 2 NG Transit Lite easier to use and push. It has our narrowest frame ever and a shorter wheelbase. Manoeuvre with less effort, even around tighter spaces. It's simply easier to get around.

*From 10.5 kg transport weight.

2



Extra comfort.

We've designed a new Y-shaped backrest for extra comfort. Smooth steering and strong brakes let you travel with confidence. A newly positioned parking brake is accessible from the seated position, helping to maintain levels of independence.

3



Folds in seconds.

The compact folding frame and semi-folding backrest lets Action 2 NG Transit Lite fit snugly into any space. Removable legrests and armrests detach in seconds. Easy to fold down and set up, always ready for your next journey.



Attendant braking system

Ease of transfer
Swing away and removable armrests and legrests.

Ideal for storage

New user brakes

New lighter front and rear wheels



What's new






















Reduced weight and frame size



Easily accessible parking brake

Technical data

For more information about this product, including the user manual, please visit our website: www.invacare.eu

Action 2 NG Transit Lite	Seat width	Seat depth	Seat angle	Seat height	Backrest angle	Backrest height	Front frame angle	
								
	405 / 430 / 455 / 480 mm	400 mm	3°	485 / 510 mm	7°	430 mm	80°	
	Action 2 NG Transit Lite	Length of footrests	Armrest height	Driving unit width	Total width, folded	Total length inc. Footrests	Total length without footrests	Total height
								
		330 / 460 mm	280 mm	Seat width +100 mm	250 mm	995 mm	760 mm	975 mm
		Action 2 NG Transit Lite	Max. Safe slope	Radius	Max. User weight	Total weight	Mass of the heaviest part	
								
7°			870 mm	125 kg	13.1 kg	10.2 kg		

Frame colours



Mocha



Discover the
full Action range



Tel: +44 61 487 70 70
hqeurope@invacare.com
www.invacare.eu.com

Invacare International GmbH
Neuhofweg 51
4147 Aesch BL
Switzerland

© 2023 Invacare International GmbH. All rights reserved. All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation. EU - 10/23



Yes, you can.®