



# GO ANYWHERE DO ANYTHING

The **Compact attract** makes daily tasks really simple. This cross-folding wheelchair is the perfect basic equipment and has an unsurpassed look. Our proven quality is reflected in this wheelchair, and it is the perfect companion for your activities. Ingenious ideas make it compact. Its precision makes it robust.



Transport size by a small width



T-armrest with soft pad and clothes guards



Timeless and easy



















Vertical adapter plate for easy adjustment of the rear wheel position



Push brake in open position

## TECHNICAL DATA COMPACT ATTRACT



-  **Frame angle**  
70°/90°
-  **Seat height**  
front: 430 – 520 mm  
rear: 400 – 490 mm  
in increments of 10 mm
-  **Max. safe slope**  
7°
-  **Total weight unloaded**  
approx. 11.3 kg
-  **Seat width**  
360 – 500 mm  
in increments of 20 mm
-  **Knee-to-heel length**  
380 – 460 mm  
in increments of 10 mm
-  **Reversing width**  
920 – 1125 mm
-  **Transport weight w/out rear wheels**  
approx. 8.7 kg
-  **Seat depth**  
360 – 500 mm  
in increments of 20 mm
-  **Total width folded**  
approx. 300 mm
-  **Total width**  
SW + 180 mm
-  **Max. user weight**  
130 kg
-  **Backrest height**  
330 – 495 mm  
in increments of 15 mm
-  **Backrest angle**  
82°/86°/90°/  
94°/98°/102°
-  **Total length**  
70°: approx. 985 mm  
90°: approx. 895 mm  
w/out legrest: 690 mm
-  **Crash Test approved**  
Fulfills the requirements of ISO 7176-19