

AQUATEC®

Sorrento



Bathtub Swivel Seat



Yes, you can.®

The *Aquatec Sorrento Bathtub Swivel Seat* is one of the outstanding modular and modern products of the *Aquatec Sorrento* range.

With great function but also with great style and design. The non-slip soft touch surfaces as well as the soft backrest provide a secure and very comfortable seat for the user. The 4 anti-slip attachment legs assure a secure positioning of the bathtub swivel seat and will cover all standard bathtubs. The all-corrosion free, tool-less assembled bathtub swivel seat is the best choice for users who need assistance in the bathroom for their personal hygiene.



Sorrento - looks great and perfect in every bathroom

Swivel function

It was never so easy to use a **Bathtub Swivel Seat**. With the release handle on both sides and improved technology of the swivel function you will be highly satisfied with the performance of the *Aquatec Sorrento Bathtub Swivel Seat*.



Non-slipping Soft Touch Surface

The *Aquatec Sorrento* is the first Bathroom collection with the innovative TPE Soft Touch Surface on the seat and the armrests. Not only the great look but also the comfortable feeling will persuade the user.

Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

	Sorrento Swivel Seat
Bathtub Swivel Seat	1532693
Measurements	
Seat width (mm)	400
Seat depth (mm)	400
Total width (mm)	710-760
Total depth (mm)	650
Total height (mm)	520
Width between armrests (mm)	450
Armrest length (mm)	280
Armrest height (mm)	195
Backrest width (mm)	500
Backrest height (mm)	380
Backrest cover height (mm)	280
Width hygiene recess (mm)	140
Depth hygiene recess (mm)	170
High load capacity	130 kg
Weight	6 kg
Colour	



Tel: +41 61 487 70 70
 hqeuropa@invacare.com
 www.invacare.eu.com

Invacare International GmbH
 Neuhofweg 51
 4147 Aesch BL
 Switzerland



Yes, you can.®