





Monteringsanvisning
 Assembly instruction
 Bedienungsanleitung
 Instructions de montage





Art no 12515

OBS! Bromsen får inte vara i bromsat läge vid montering av handtaget.
 NB! The brake must not be in active position when assembling the handle.
 ACHTUNG! Bei Montage des Handgriffes muss die Bremse gelöst sein.
 Attention ! Le frein ne doit pas être en position de freinage lors du montage de la poignée.

Demontage Disassembly Demontage Démontage



Montering Assembly Montage Montage



T20
 Torx-nyckel
 Torx-key
 Torx-schlüssel
 Clef mâle pour
 six-pans.

